


FINGER PRICK SAMPLE COLLECTION:


- Wipe the complete finger tip with the alcohol swab as shown in Fig. (d). Wait until the finger has completely dried (minimum 30 seconds).
- Take the lancet and prick the side of the pulp (ball of the finger) with the lancet, perpendicular to the lines of the finger print as shown in Fig. (e).
- Make sure a well formed drop of blood is present on the tip of the finger.
- Take the sample dropper and collect 4µl of blood by dipping the tip of the sample dropper into the blood drop as shown in Fig. (f) and immediately place the tip of the sample dropper in the sample well "S". (Care should be taken that the blood sample does not clot & the transfer to the sample pad is immediate).
- Add whole blood sample using the sample dropper/ micropipette onto the sample pad in the sample well 'S'. fig. (g)




- NOTE :** Press the tip of the dropper onto the sample pad in the sample well "S" to ensure that the complete volume of whole blood has been transferred to the strip.
- Add 3 drops of the Assay Buffer in the buffer well 'B'. fig. (h). Ensure FREE FALLING OF DROPS on the membrane, holding the vial/dropper vertically for proper volume. Screw cap the vial after use.
 - Allow the reaction to occur for 20 minutes.
 - Read the results at 20 minutes. Do not read the result after 20 minutes. Reading beyond prescribed time may give false results.
 - Discard the ADVANTAGE MAL CARD immediately after reading results at 20 minutes as it is potentially infectious.

INTERPRETATION OF THE RESULTS

REACTIVE




As shown in fig. (i), appearance of three purplish pink coloured lines one each in Pf. region (F), Pan region (P) & Control region (C) indicates that the sample is reactive for P. falciparum or mixed infection of Pf and P. vivax or (P. malariae / P. ovale).







As shown in fig. (j) appearance of two purplish pink coloured line one each at P & C region only indicates that the sample is reactive for P. vivax / P. malariae / P. ovale only. A difference of intensity in colour may be observed between both the test lines ('P' & 'F') and between the test lines and control line from faint to strong intensity depending on the concentration of pLDH antigen in the sample. In case of faint test line ('P' and/or 'F') repeat the test using fresh card.

NON-REACTIVE



As shown in fig. (k), appearance of only one purplish pink coloured line at Control (C) region indicates that the sample is non-reactive for P. falciparum and other Plasmodium Species (P. vivax / P. malariae / P. ovale).

INVALID



The test is invalid, if or no control line appears after the completion of test either with clear background or with complete pinkish/ purplish background fig. (l, m, n & o). The test should be repeated using a new card.

LIMITATIONS AND INTERFERENCES

- The test procedure, precautions and interpretation of results for this test must be strictly followed.
- As with all diagnostic tests, the test result must always be correlated with clinical finding.

- Though the test is accurate in detecting pLDH specific to P. falciparum and other Plasmodium species (P. vivax / P. malariae / P. ovale), a low incidence of false results can occur. All reactive test results are to be interpreted within the epidemiological, clinical and therapeutic context. When it seems indicated, the parasitological techniques of reference should be considered (microscopic examination of the thick smear and thin blood films).
- Any modification to the above procedure and / or use of other reagents will lead to erratic results.
- In case of faint test line, repeat the test using fresh card. If faint test line persists, confirm the test result with ELISA/Cia/ microscopy and co-relate clinically.
- In most of the cases, after successful anti-malarial therapy the Pan band will turn negative. However, depending on the medication used, the clearance of parasite may take longer and the reaction of the test may remain positive. In such cases the test should be repeated after 5-10 days of start of treatment.
- In case only "F" test line appears and "P" line is absent on the card for a sample, then confirm the result with microscopy. If the sample is negative by microscopy then interpret the result as Malaria Negative. This could be due to cross reactivity with Rheumatoid factor and heterophilic antibodies.
- The possibility of resistant strain of malaria should always be considered if the reaction of the test remains positive with the same intensity after 5-10 days post treatment.
- Patient with rheumatoid factors, anti-nuclear antibody or dengue may give false positive results.

PERFORMANCE CHARACTERISTICS OF ADVANTAGE MAL CARD

(i) WHO Evaluation:

The ADVANTAGE MAL CARD test kit has been evaluated by WHO, Geneva using a panel of wild & cultured malaria positive samples and the results obtained are as follows :

| | Panel Detection Score* | Specificity* |
|------------------|---------------------------|--------------|
| 200 Parasites | 94.3% (Pv) | — |
| 2000 Parasites | 94.0% (Pf) and 97.1% (Pv) | — |
| Negative Samples | — | 99.6% |

*Reference: Malaria Rapid Diagnostic Test Performance: Results of WHO product testing of malaria RDTs: Round 5 (2013), Page:38.

Note: The above information is provided for the scientific community, It is not for commercial or promotional purpose.

(ii) In-house Evaluation:

Analytical Sensitivity: The test can detect parasitemia levels of ≥ 100 parasites per μl of blood for P. falciparum (pLDH) and ≥ 200 parasites per μl of blood for P. vivax (pLDH).

The ADVANTAGE MAL CARD has been evaluated in-house with malaria positive and negative clinical whole blood samples and compared with microscopic examination. The evaluation also included cross-reacting samples; Dengue, Rheumatoid factor, Leptospira, HIV, HCV, HBV, M. tuberculosis, Syphilis, Brucella, Scrub typhus positive samples. The results obtained are as follows:

| Sample | Total no. of samples tested | ADV. MAL CARD | | Sensitivity (%) | Specificity (%) |
|------------------------|-----------------------------|---------------|----------|-----------------|-----------------|
| | | Positive | Negative | | |
| Malaria Negative | 2010 | 1 | 2009 | - | 99.95 |
| Cross-reacting sample | 59 | 0 | 59 | - | 100 |
| P. falciparum Positive | 58 | 58 | 0 | 100 | - |
| P. vivax Positive | 105 | 105 | 0 | - | 100 |

Precision: Within-run and between-run precisions have been determined by testing 10 replicates of 8 specimens: four negative and four positive; two weak and two moderate positive. The C.V. (%) of negative, moderate positive and weak positive samples were within 10% of the time.

LIMITED EXPRESSED WARRANTY DISCLAIMER

The manufacturer limits the warranty to the test kit, as much as that the test kit will function as an in vitro diagnostic assay within the limitations and specifications as described in the product instruction-manual, when used strictly in accordance with the instructions contained therein. The manufacturer disclaims any warranty expressed or implied including such expressed or implied warranty with respect to merchantability, fitness for use or implied utility for any purpose. The manufacturer's liability is limited to either replacement of the product or refund of the purchase price of the product and in no case liable to claim of any kind for an amount greater than the purchase price of the goods in respect of which damages are likely to be claimed. The manufacturer shall not be liable to the purchaser or third parties for any injury, damage or economic loss, howsoever caused by the product in the use or in the application thereof.

BIBLIOGRAPHY OF SUGGESTED READING

- Jamshaid Iqbal; etal. Journal of Clinical Microbiology, 37(1), 1999, 3644-3646.
- World Health Organization 2000. New perspectives malaria diagnosis. World Health Organization, Geneva, Switzerland.
- Piper R, etal. Am J. Trop. Med. Hyg., 60(1), 1999, 109-118.
- Moody A, Clinical Microbiological Review, 15(1), 2002, 66-78.

in vitro diagnostic reagent, not for medicinal use

J. MITRA & CO. PVT. LTD.

A 180-181, Okhla Indl. Area, Phase-1, New Delhi-110 020, INDIA

Ph: +91-11-47130300, 47130500

e-mail: jmitra@jmitra.co.in Internet: www.jmitra.co.in