


TEST PROCEDURE

RT
20-30°C

1. Bring the complete kit and specimen to be tested to room temperature prior to testing.
2. Remove the test cartridge from the foil pouch prior to use and place it on a flat and dry surface. The test should be performed immediately after removing the test card from the foil pouch.
3. Label the test cartridge with patient's name or identification number. **Do not write on QR code.**
4. Add 50µl Serum/ Plasma sample using micropipette onto the sample pad in the sample well 'S'. **Care should be taken to avoid any spillage on the QR Code.**



NOTE : Make sure that the sample from the micropipette has been completely transferred to the sample pad.

5. Add 2 drops of the Assay Buffer in the buffer well 'B'.
6. Allow the reaction to occur for 30 minutes. In the meantime enter the patient's details in the iQuant analyzer testing window and select the NS1 test from the pop down menu in the testing window of the iQuant analyzer.
7. Insert the test cartridge into the i-Quant Analyzer with arrow (←) marked side on the top of cartridge facing towards the analyzer and press RUN icon. Note down the value displayed on the screen of i-Quant Analyzer & interpret the result as mentioned below.
8. Discard the Dengue NS1 Ag Quanti Card immediately after reading results at 30 minutes considering it to be potentially infectious.

Important Note: Do not read results after 30 minutes.

INTERPRETATION OF RESULTS

The iQuant Analyzer will display results as Reactive, Equivocal or Non-Reactive as follows:

NS1 Antigen Reactive: >1.1 U and above: interpret the result as Dengue NS1 Antigen reactive.

Equivocal: >0.9 to ≤1.1 U: Interpret the result as equivocal. Repeat the test after centrifuging the sample at 5000 rpm for 20 minutes. Even after repeating the test, if result comes equivocal, further test the sample with alternative method or collect another sample.

Non-Reactive: Below ≤0.9 U: Interpret the result as Dengue NS1 Antigen non-reactive.

LIMITATIONS AND INTERFERENCES

1. The test is for in vitro diagnostic use only.
2. This test detects the presence of Dengue NS1 antigen in the specimen and should not be used as the sole criteria for the diagnosis of Dengue virus infection.
3. Serological cross-reactivity across the Flavivirus group (Dengue virus, St. Louis encephalitis, Japanese encephalitis, West Nile and yellow fever virus) is common.
4. As with all diagnostic tests, all results must be correlated with other clinical findings. If the test result is negative and clinical symptoms persist, additional follow-up testing using other clinical methods is recommended. A negative result at any time does not preclude the possibility of an early infection of Dengue virus.
5. This is only a screening test. Therefore, isolation of virus, antigen detection in fixed tissues, RT-PCR and more specific alternative diagnosis method must be used in order to obtain a confirmation of dengue virus infection.

PERFORMANCE CHARACTERISTICS OF DENGUE NS1 Ag QUANTI CARD TEST

The kit has been evaluated in-house with the known panel of fresh as well as frozen Dengue NS1 antigen positive and Negative samples. The performance of the test kit was evaluated and compared with the license commercially available ELISA test kit. The samples

included cross-reacting samples; Epstein-Barr virus, Malaria, Rheumatoid factor, Leptospirosis, Japanese encephalitis, yellow fever and West Nile viruses. Following is the in-house evaluation.

| Sample Type | No. of Samples tested | Result of licensed test | Dengue NS1 Ag Quanti Card |
|-----------------------------|-----------------------|-------------------------|---------------------------|
| Dengue NS1 Antigen Negative | 2070 | 2070 | 2069 |
| Dengue Antigen Positive | 255 | 255 | 255 |

Sensitivity : **100%**

Specificity : **99.95%**

Precision:

Intra assay precision (Reproducibility)

Within run (Intra assay) precision have been determined by testing 10 replicates of two negative and four dengue NS1 antigen positive samples; 3 weak and 1 medium. The C.V. (%) of all the samples were within 10% of the unit.

Inter assay precision (Reproducibility)

Between run (Inter assay) precision have been determined by testing 10 replicates of two negative and four dengue NS1 antigen positive samples; 3 weak and 1 medium in 10 sequential days. The C.V. (%) of all the samples were within 10% of the unit.

LIMITED EXPRESSED WARRANTY DISCLAIMER

The manufacturer limits the warranty to the test kit, as much as that the test kit will function as an in vitro diagnostic assay within the limitations and specifications as described in the product instruction-manual, when used strictly in accordance with the instructions contained therein. The manufacturer disclaims any warranty expressed or implied including such expressed or implied warranty with respect to merchantability, fitness for use or implied utility for any purpose. The manufacturer's liability is limited to either replacement of the product or refund of the purchase price of the product and in no case liable to claim of any kind for an amount greater than the purchase price of the goods in respect of which damages are likely to be claimed. The manufacturer shall not be liable to the purchaser or third parties for any injury, damage or economic loss, howsoever caused by the product in the use or in the application thereof.

TROUBLE SHOOTING

| PROBLEM | POSSIBLE CAUSE | SOLUTION |
|-------------------|---|--|
| 1. False negative | a) Low sample volume along with less volume of Assay Buffer used. | Repeat the test with proper volume of samples and/or Assay Buffer. |
| 2. False positive | a) Less volume of Assay Buffer added/used. | Repeat the test with proper volume of Assay Buffer |

BIBLIOGRAPHY OF SUGGESTED READING

1. Dengue haemorrhagic fever: Diagnosis, Treatment, Prevention and Control. WHO 2009 edition.
2. Evaluation of diagnostic test: Dengue, Rosanna W. Peeling, Harrey Artsob etal. (2010). Nature reviews.
3. Use of dengue NS1 antigen for early diagnosis of dengue virus infection. Kassim FM; Izzate MN, etal. (2011) Southeast Asian J. Trop. Med. Public Health. May; 42(3); 562-9.
4. Dengue NS1 antigen detection: A useful tool in early diagnosis of dengue virus infection. S. Datta, C Wattal (2010). Indian Journal of Medical Microbiology, vol 28, No. 2; 107-110.
5. Dengue: guidelines for diagnosis, treatment, prevention and control –New Edition WHO/HTM/NTD/DEN/2009.1

For in-vitro diagnostic use only, not for medicinal use

J. MITRA & CO. PVT. LTD.

A 180-181, Okhla Indl. Area, Phase-1, New Delhi-110 020, INDIA

Ph: +91-11-47130300, 47130500, 26818971-73

e-mail: jmitra@jmitra.co.in Internet: www.jmitra.co.in