

DENGUE NS1 Ag MICROLISA

Microwell ELISA Test for the Detection of Dengue NS1 Antigen in Human Serum/Plasma

1. INTRODUCTION

Dengue virus is a flavivirus found largely in areas of the tropic and sub-tropics. There are four distinct but antigenically related serotypes of dengue viruses, and transmission is by mosquito, principally *Aedes aegypti* and *Aedes albopictus*.

The mosquito-borne dengue viruses (serotype 1-4) cause dengue fever, a severe flu-like illness. The disease is prevalent in third world tropical regions and spreading to sub-tropical developed countries - including the United States. WHO estimates that 50-80 million cases of dengue fever occur worldwide each year, including a potentially deadly form of the disease called dengue haemorrhagic fever (DHF) and dengue shock syndrome (DSS). Primary infection with dengue virus results in a self-limiting disease characterized by mild to high fever lasting 3 to 7 days, severe headache with pain behind the eyes, muscle and joint pain, rash and vomiting. Secondary infection is the more common form of the disease in many parts of Southeast Asia and South America. This form of the disease is more serious and can result in DHF and DSS. The major clinical symptoms can include high fever, haemorrhagic events, and circulatory failure, and the fatality rate can be as high as 40%. Early diagnosis of DSS is particularly important, as patients may die within 12 to 24 hours if appropriate treatment is not administered.

Primary dengue virus infection is characterized by elevations in specific NS1 antigen levels 0 to 9 days after the onset of symptoms; this generally persists upto 15 days. Earlier diagnosis of Dengue reduces risk of complication such as DHF or DSS, especially in countries where dengue is endemic.

2. INTENDED USE

















DENGUE NS1 Ag MICROLISA is designed for *in vitro* qualitative detection of Dengue NS1 antigen in human serum or plasma and is used as a screening test for testing of collected blood samples suspected for DENGUE. The kit detects all four subtypes; DEN1, DEN2, DEN3 & DEN4 of Dengue Virus.

3. PRINCIPLE

DENGUE NS1 Ag MICROLISA is a solid phase enzyme linked immunosorbent assay (ELISA) based on the "Direct Sandwich" principle. The microwells are coated with Anti-dengue NS1 antibodies with high reactivity for Dengue NS1 Ag. The samples are added in the wells followed by addition of enzyme conjugate (monoclonal anti-dengue NS1 antibodies linked to Horseradish Peroxidase (HRPO)). A sandwich complex is formed in the well wherein dengue NS1 (from serum sample) is "trapped" or "sandwiched" between the antibody and antibody HRPO conjugate. Unbound conjugate is then washed off with wash buffer. The amount of bound peroxidase is proportional to the concentration of dengue NS1 antigen present in the sample. Upon addition of the substrate buffer and chromogen, a blue colour develops. The intensity of developed blue colour is proportional to the concentration of dengue NS1 antigen in sample. To limit the enzyme-substrate reaction, stop solution is added and a yellow colour develops which is finally read at 450nm spectrophotometrically.

4. DESCRIPTION OF SYMBOLS USED














The following are graphical symbols used in or found on J. Mitra diagnostic products and packing. These symbols are the most common ones appearing on medical devices and their packing. They are explained in more detail in the European Standard EN ISO 15223-1:2021.

 Manufactured By	 <i>In vitro</i> diagnostic medical device
 No. of tests	 See Instruction for use
 Lot Number Batch Number	 Temperature Limitation
 Manufacturing Date	 Caution See instruction for use
 Expiry Date	 Catalogue Number
 Do not use if package is damaged	 Keep away from sunlight
 Contains biological Material of Human Origin	 Contains biological Material of Animal Origin
 Country of Manufacture	 Keep Dry

5. KIT PRESENTATION

- 96 Test Pack

6. KIT & ITS COMPONENTS

COMPONENT	DESCRIPTION
 Microwells	Microwells coated with anti-Dengue NS1 antibodies packed in a sealed pouch with desiccant.
 Diluent	Buffer containing protein stabilizers & antimicrobial agents as preservative and to be used for Sample & Conjugate dilution.
 Enzyme Conjugate Concentrate (50X)	Containing Monoclonal Anti-Dengue NS1 linked to horseradish peroxidase with protein stabilizers.
 Wash Buffer Concentrate (25X)	Concentrated Phosphate buffer with surfactant.
 TMB Substrate	TMB, to be diluted with TMB Diluent before use.
 TMB Diluent	Buffer solution containing H ₂ O ₂ with preservative.
 Control - 	Normal human serum negative for Dengue NS1 antigen with preservative.
 Control + 	Recombinant Dengue NS1 antigen, with preservative.
 Calibrator	Recombinant Dengue NS1 antigen, with preservative.
 Stop Solution	Ready to use, 1N H ₂ SO ₄ .
 Plate Sealers	Adhesive sheets to cover the microwells during incubation.

7. STORAGE AND STABILITY

Store the kit & its components at 2-8°C. Expiry date on the kit indicates the date beyond which kit should not be used.


8. ADDITIONAL MATERIAL AND INSTRUMENTS REQUIRED

- Micropipettes and microtips
- Timer
- Elisa reader
- Elisa washer
- Distilled or deionized water
- Incubator 37°C
- Graduated Cylinders, for reagent dilution
- Disinfectant Solution
- Paper towels or absorbent tissue
- Disposable gloves

9. SPECIMEN COLLECTION & HANDLING

1. Only human serum or plasma samples should be used for the test. While preparing serum samples, remove the serum from the clot as soon as possible to avoid hemolysis. Fresh serum/plasma samples are preferred.
2. Specimens should be free of microbial contamination and may be stored at 2-8°C for one week, or frozen at -20°C or lower. Avoid repeated freezing and thawing.
3. Do not use heat inactivated samples as their use may give false results. Hemolyzed and Icteric hyperlipemic samples may give erroneous results.
4. Do not use Sodium Azide as preservative because it inactivates Horseradish peroxidase.

10. WARNING & PRECAUTION

 **CAUTION:** THIS KIT CONTAINS MATERIALS OF HUMAN ORIGIN. NO TEST METHOD CAN OFFER COMPLETE ASSURANCE THAT HUMAN BLOOD PRODUCTS WILL NOT TRANSMIT INFECTION. NEGATIVE CONTROL, POSITIVE CONTROL & ALL THE SAMPLES TO BE TESTED SHOULD BE HANDLED AS THOUGH CAPABLE OF TRANSMITTING INFECTION.

1. The use of disposable gloves and proper biohazardous clothing is STRONGLY RECOMMENDED while running the test.
2. In case there is a cut or wound in hand, DO NOT PERFORM THE TEST.
3. Do not smoke, drink or eat in areas where specimens or kit reagents are being handled.
4. Tests are for *in vitro* diagnostic use only and should be run by competent person only.
5. Do not pipette by mouth.
6. All materials used in the assay and samples should be decontaminated in 5% sodium hypochlorite solution for 30-60 min. before disposal or by autoclaving at 121°C at 15psi for 60 min. Do not autoclave materials or solution containing sodium hypochlorite. They should be disposed off in accordance with established safety procedures.
7. Wash hands thoroughly with soap or any suitable detergent, after the use of the kit. Consult a physician immediately in case of accident or contact with eyes, in the event that contaminated material are ingested or come in contact with skin puncture or wounds.

8. Spills should be decontaminated promptly with Sodium Hypochlorite or any other suitable disinfectant.
9. Stop solution contains sulfuric acid. If sulfuric acid comes in contact with the skin, wash thoroughly with water. In case of contact with eyes, flush with excess of water.
10. ELISA Reader & micropipettes used in testing should be calibrated at regular interval to ensure accurate results.

11. SPECIMEN PROCESSING

(A) FROZEN SAMPLE

Dengue NS1 Ag Microlisa test is best used with fresh samples that have not been frozen and thawed. However most frozen samples will perform well if the procedure suggested below is followed.

Allow the frozen sample to thaw in a vertical position in the rack. Do not shake the sample. This allows particles to settle to the bottom. If a centrifuge is available, the sample should be centrifuged. (10,000 rpm for 15 min.)

(B) TRANSPORTATION

If the specimen is to be transported, it should be packed in compliance with the current Government regulations regarding transport of aetiologic agents.

12. PRECAUTIONS FOR USE

Optimal assay performance requires strict adherence to the assay procedure described in the manual.

1. Do not use kit components beyond the expiration date, which is printed on the kit.
2. Avoid microbial contamination of reagents. The use of sterile disposable tips is recommended while removing aliquots from reagent bottles.
3. Prepare working substrate solution just 10 minutes prior to adding in the wells.
4. If blue colour or white particles appears in working substrate solution then do not use it. Take fresh containers and tips and prepare it again.
5. Use separate tips for TMB substrate and TMB diluent.
6. Mark the test specimen with patient's name or identification number. Improper identification may lead to wrong result reporting.
7. Do not allow microwells to dry once the assay has started.
8. Ensure that the microwell strips are levelled in the strip holder. Before reading, wipe the bottom of the microwell strips carefully with soft, absorbent tissue to remove any moisture.
9. If available, a microwell reader which contains a reference filter with settings at 620 or 630 nm should be used. Use of a reference filter minimises interference due to microwells that are opaque, scratched or irregular. However, if a reference filter is unavailable, the absorbance may be read at 450 nm without a reference filter.
10. Distilled or deionised water must be used for wash buffer preparation.
11. Bring all the reagents to room temperature (20-30°C) before use.
12. Do not combine reagents from different batches, as they are optimized for individual batch to give best results.
13. Due to interchange of caps the reagents may get contaminated. Care should be taken while handling the reagents to avoid any sort of contamination.
14. Run negative control, positive control and calibrator in each assay.
15. Use freshly collected, clean serum/ plasma samples for assay. Try to avoid Haemolyzed/ turbid/ lipemic serum or plasma samples.
16. Use a separate tip for each sample and then discard it as biohazardous waste.
17. Thorough washing of the wells is critical to the performance of the assay.
18. Avoid strong light exposure during the assay.

13. PRELIMINARY PREPARATIONS

- Pre-warm the incubator to + 37°C.
 - Bring foil pack to room temperature (20-30°C) before opening to prevent condensation on the microwell strips.
- a. Break-off the required number of strips needed for the assay and place in the strip holder. Take the strip holder with the required number of strips, taking into account that, one negative control, one positive control and three calibrator should be included in the run while opening the fresh kit. However for one or two strips one negative control, one positive control and two calibrator should be included in each subsequent runs.
 - b. **Unused wells should be stored at 2-8°C with desiccant in an aluminium pouch with clamp & rod. Microwells are stable for 30 days at 2-8°C from the date of opening of sealed pouch, when stored with desiccant along with clamp & rod.**

Caution: Handle microwell strip with care. Do not touch the bottom exterior surface of the wells.

● Preparation of Working Wash Buffer:

- a) Check the buffer concentrate for the presence of salt crystals. If crystals are present in the solution, resolubilize by warming at 37°C until all crystals dissolve.
- b) Prepare at least 25ml. (1ml. concentrated buffer with 24 ml. water) of buffer for each strip used. Mix well before use.
- c) Mix 20 ml. of 25X wash buffer concentrate with 480ml. of distilled or deionized water. Wash buffer is stable for 2 months when stored at 2-8°C.

● Preparation of Working Conjugate:

Dilute conjugate concentrate 1:50 in Diluent. **Do not store working conjugate.** Prepare a fresh dilution for each assay in a clean glass vessel. Determine the quantity of working conjugate solution to be prepared from the table below. Mix solution thoroughly before use.

No. of Strips	1	2	3	4	5	6	7	8	9	10	11	12
No. of Wells	8	16	24	32	40	48	56	64	72	80	88	96
Enzyme Conjugate Concentrate (µl)	20	40	60	80	100	120	140	160	180	200	220	240
Diluent (ml)	1.0	2.0	3.0	4.0	5.0	6.0	7.0	8.0	9.0	10.0	11.0	12.0

● Preparation of working substrate solution :

Mix TMB substrate and TMB Diluent in 1:1 ratio to prepare working substrate.

No. of Strips	1	2	3	4	5	6	7	8	9	10	11	12
No. of Wells	8	16	24	32	40	48	56	64	72	80	88	96
TMB Substrate (ml)	0.8	1.6	2.4	3.2	4.0	4.8	5.6	6.4	7.2	8.0	8.8	9.6
TMB Diluent (ml)	0.8	1.6	2.4	3.2	4.0	4.8	5.6	6.4	7.2	8.0	8.8	9.6

Do not store working substrate. Prepare a fresh dilution for each assay in a clean plastic/glass vessel. Determine the quantity of working substrate solution to be prepared from table. Mix solution thoroughly before use. Discard unused solution. A deep blue color present in the substrate solution indicates that the solution has been contaminated and must be discarded.

15. PROCEDURAL NOTES:

1. Material should not be used after the expiry date shown on the labels. Components and test specimen should be at room temperature (20-30°C) before testing begins. Return the reagents to 2-8°C after use.
2. Mark the test specimen with patient's name or identification number. Improper identification may lead to wrong result reporting.
3. All reagents must be mixed well before use.
4. To avoid contamination, do not touch the top or bottom of strips or edge of wells.
5. All pipetting steps should be performed with utmost care and accuracy. Cross contamination between reagents and samples will invalidate results.
6. Prevent evaporation during sample incubation by covering the strips with sealer; remove sealer before washing.
7. Routine maintenance of wash system is strongly recommended to prevent carry over from highly reactive specimens to non reactive specimens.











16. TEST PROCEDURE

Once the assay has started, complete the procedure without interruption. All the reagents should be dispensed in the centre of the well and the tip of the pipette should not touch the wall of the microwell.

Fit the strip holder with the required number of Anti-Dengue NS1 antibody coated strips. The sequence of the procedure must be carefully followed. Arrange the assay control wells in a horizontal or vertical configuration. Configuration is dependent upon reader software.

1. Add 50 µl Diluent in all the wells.
2. Add 50 µl Negative Control in A-1well.
3. Add 50 µl Calibrator in B-1, C-1 & D-1 well.
4. Add 50 µl Positive Control in E-1 well.
5. Add 50 µl sample in F-1 well onwards.
6. Add 100 µl of working Conjugate Solution in each well.
7. Ensure thorough mixing of controls, samples to be tested & working conjugate to get reproducible results.
8. Apply cover seal.
9. Incubate at 37°C ± 1°C for 90 min. ± 1min.
10. While the samples and working Conjugate are incubating, prepare working Wash Solution as specified in preparation of reagents.
11. Take out the plate from the incubator after the incubation time is over and, wash the wells 6 times with working Wash Solution.
12. Add 150 µl of working substrate solution in each well.

13. Incubate at room temperature (20-30°C) for 30 ± 2 minutes in dark.
14. Add 100 µl of stop solution and leave at room temperature (20-30°C) for 1 minute.
15. Read absorbance at 450 nm within 5 minutes of adding stop solution in an ELISA READER. (Bichromatic absorbance measurement with a reference wavelength 600-650 nm is recommended when available).

SUMMARY OF PROCEDURE		
Add controls, calibrator and samples		50 µl Diluent 50 µl controls, calibrator and samples
Prepare working conjugate		No of Strips 1 2 3 4 5 6 7 8 9 10 11 12 Enz. conc. 20 40 60 80 100 120 140 160 180 200 220 240 (µl) Diluent (ml.) 1 2 3 4 5 6 7 8 9 10 11 12
Add Conjugate		100 µl
Cover the plate & incubate		90 minutes at 37°C
Wash		6 Cycles
Prepare Chromogenic Substrate		No of Strips 1 2 3 4 5 6 7 8 9 10 11 12 TMB 0.8 1.6 2.4 3.2 4.0 4.8 5.6 6.4 7.2 8.0 8.8 9.6 (ml) Substrate 0.8 1.6 2.4 3.2 4.0 4.8 5.6 6.4 7.2 8.0 8.8 9.6 (ml.) Diluent (ml.)
Add Substrate		150 µl
Incubate in dark		30 ± 2 minutes at Room Temp.
Add Stop Solution		100 µl
Read Results		450 nm / 630 nm

17. CALCULATION OF RESULTS

TEST VALIDITY:

Ensure the following is within specified acceptance criteria

- i) NC O.D. must be < 0.3. If it is not so, the run is invalid and must be repeated.
- ii) PC O.D. must be > 1.0. If it is not so, the run is invalid and must be repeated.
- iii) Mean Calibrator O.D. must be ≥ 0.35. If it is not so, the run is invalid and must be repeated.
- iv) Cut off value must be ≥ 1.5 x NC O.D. If it is not so, the run is invalid and must be repeated.
- v) Ratio of PC O.D. / cut off must be > 1.1. If it is not so, the run is invalid and must be repeated.

Imp. Note: The calibration factor detail is batch specific and stamped on back page of Instruction manual.

- a. Cut off value = mean O.D. of calibrator x calibration factor
- b. Calculation of sample O.D. ratio : Calculate sample O.D. ratio as follows:

$$\text{Sample O.D. ratio} = \frac{\text{Sample O.D.}}{\text{Cut off Value}}$$

- c. Calculation of Dengue NS1 Ag units : Calculate by multiplying the sample O.D. ratio by 10.
Dengue NS1 Ag units = sample O.D. ratio x 10.
e.g.: Mean O.D. of calibrator = 0.75
Calibration factor = 0.7
Cut off value = 0.75 x 0.7 = 0.525
e.g.: sample absorbance (O.D.) = 0.925
Cut off value = 0.525
Sample O.D. ratio = 0.925 / 0.525 = 1.761
Dengue NS1 Ag units = 1.761 x 10 = 17.61

18. INTERPRETATION OF RESULTS

- a. If the Dengue NS1 Ag Units is < 9 then interpret the sample as Negative for Dengue NS1 Antigen.
- b. If the Dengue NS1 Ag Units is between 9 - 11 then interpret the sample as Equivocal for Dengue NS1 Antigen.
- c. If the Dengue NS1 Ag Units is > 11 then interpret the sample as Positive for Dengue NS1 Antigen.

19. PERFORMANCE CHARACTERISTICS

External Evaluation:

Sensitivity with RT-PCR confirmed positive samples

The sensitivity of Dengue NS1 Ag Microlisa kit has been evaluated on 58 retrospective serum samples from patients with dengue infection confirmed by Dengue RT-PCR. The samples were from the four different dengue serotypes (DEN 1, DEN 2, DEN 3 & DEN 4).

The sensitivity of the Dengue NS1 Ag Microlisa kit when compared with Dengue RT- PCR is as follows:

Serotype	Dengue RT PCR Positive Samples	Result of Dengue NS1 Ag Microlisa	Sensitivity of Dengue NS1 Ag Microlisa
Den 1	25	25	100%
Den 2	8	7	87.5%
Den 3	22	22	100%
Den 4	3	3	100%

The overall sensitivity of Dengue NS1 Ag Microlisa kit on the above panel was found to be 98.28%.

In-house Evaluation:

The kit has also been evaluated in-house with the known panel of fresh as well as frozen Dengue NS1 antigen positive of all 4 dengue serotypes and Negative samples. The samples included cross-reacting samples; Epstein-Barr virus, Malaria, Rheumatoid factor, Leptospirosis, Japanese encephalitis, yellow fever and West Nile viruses. Following are the result of evaluation:

No. of Samples	Status	Dengue NS1 Ag Microlisa	
		Positive	Negative
200	Dengue NS1 Ag Positive	199	1
2050	Dengue Negative	0	2050

Sensitivity : 99.5%

Specificity : 100%

Precision : Within-run and between-run precisions have been determined by testing 10 replicates of five samples: 1 negative and 4 dengue NS1 Ag positive; 2 weak positive, 1 medium positive & 1 strong positive. The C.V.(%) of negative, weak positive, medium positive & strong positive values were within 10%.

20. LIMITATION OF THE TEST

1. The test should be used for detection of NS1 Ag in serum or plasma only and not in other body fluids.
2. **This is only a screening test** and will only indicate the presence or absence of Dengue NS1 antigen in the specimen. All reactive samples should be confirmed by confirmatory test. Therefore for a definitive diagnosis, the patients clinical history, symptomatology as well as serological data should be considered. The results should be reported only after complying with the above procedure.
3. False positive results can be obtained due to cross reaction with Murray Valley and encephalitis, Japanese encephalitis, yellow fever and West Nile viruses. This occurs in less than 1% of the sample tested.

21. LIMITED EXPRESSED WARRANTY DISCLAIMER

The manufacturer limits the warranty to the test kit, as much as that the test kit will function as an *in vitro* diagnostic assay within the limitations and specifications as described in the product instruction-manual, when used strictly in accordance with the instructions contained therein. The manufacturer disclaims any warranty expressed or implied including such expressed or implied warranty with respect to merchantability, fitness for use or implied utility for any purpose. The manufacturer's liability is limited to either replacement of the product or refund of the purchase price of the product and in no case liable to for claim of any kind for an amount greater than the purchase price of the goods in respect of which damages are likely to be claimed.

The manufacturer shall not be liable to the purchaser or third parties for any injury, damage or economic loss, howsoever caused by the product in the use or in the application thereof.

22. REFERENCES

1. Gubler DJ, Trent DW: Emergence of epidemic dengue/dengue hemorrhagic fever as a public health problem in the Americas. *Infect Agents Dis* 2:383-393, 1993.
2. Enzyme-linked immunoassay for dengue virus NS1 antigen in serum and filter paper blood. Tran TN, de Vries PJ, Hoang LP, Phan GT Le HQ, Tran BQ, Vo CM, Nguyen NV, Kager PA, Nagelkerke N, Groen J. *BMC infect Dis*. 2005 Jan 25; 6:13.
3. Librarty, D. H., Young, P.R., Pickering D., Endy, T. P., Kalayanarooj, S., Green, S, Vaugh, D. W., Nisalak, A., Ennis, F. A. and Rothman, A. L. (2002). High circulating levels of the dengue virus non-structural protein NS1 early in dengue illness correlate with the development of dengue haemorrhagic fever. *J. Infect. Dis*. 186:1165-1165
4. Andrew KI falconar and Paul R. Young. (1990). Immuno affinity purification of native dimer forms of flavivirus non structural glyco protein NS1.
5. Young, P. R., Hilditch, P. A., Bletchly, C., alloran, W. (2000). An antigen capture enzyme-linked immunosorbent assay specific to dengue virus protein NS1 in the sera of infected patients. *J. Clin. Microbiol*. 38:1053-1057.

23. TROUBLE SHOOTING CHART

PROBLEM	POSSIBLE CAUSE	SOLUTION
1. No colour developed at the end of assay	<p>a) Any one reagent has been added in wrong sequence.</p> <p>b) Inactivated enzyme conjugate conc. due to improper storage</p> <p>c) Microplate inactivated, due to improper storage</p> <p>d) Inactivated substrate, improper storage or preparation procedure and repeat assay.</p> <p>e) Omission of any step in test procedure</p> <p>f) Incorrect (low) incubator temperature, timing or pipetting</p> <p>g) Improper preparation of wash buffer, error of dilution, improper mixing of reagents.</p> <p>h) Kit deterioration</p>	<p>Follow the procedure meticulously & repeat assay.</p> <p>Check storage of enzyme conjugate conc. at 2-8°C and it should be free of any contamination.</p> <p>Keep unused strips in aluminium poly pouch with the dessicant pouch inside and properly closed with clamp & rod.</p> <p>Use freshly prepared substrate solution and it should be free of any contamination. Recheck</p> <p>Follow the procedure meticulously & repeat assay.</p> <p>Check incubator temperature, procedure & repeat assay</p> <p>Check procedure & repeat assay</p> <p>Check storage of the kit and should be stored at 2-8°C.</p>
2. High O.D. value of Negative control	<p>a) Plate not stopped after 30 minutes of adding stop solution</p> <p>b) Same microtip used for calibrator, positive and negative controls</p> <p>c) Nonspecific attachment/ binding of other reagent</p>	<p>Follow the procedure meticulously & repeat assay.</p> <p>Change micropipette tips while addition of calibrator/ negative/ positive control</p> <p>If plates get scratches/ aberrations during washing, non specific proteins may bind while addition of next step.</p>
3. Too much colour in all wells of the plate (high background)	<p>a) Contaminated substrate</p> <p>b) Contaminated washing solution (1X). Poor quality of water used for diluting wash buffer conc.</p> <p>c) Over incubation of substrate and delay in addition of stop solution.</p> <p>d) Insufficient washing.</p> <p>i) Washing not consistent</p> <p>ii) Filling volume not sufficient.</p> <p>iii) Insufficient no. of wash cycles.</p> <p>iv) Contaminated wash device</p> <p>e) Use of wash solution from other manufacturer.</p> <p>f) Working substrate not protected from light</p>	<p>Check substrate (TMB Diluent) it should be colourless. If blue in colour then discard and use clean disposable container.</p> <p>Check the container and quality of water used for dilution. Use of glass distilled water is preferred.</p> <p>Follow the procedure meticulously.</p> <p>Check wash device, fill the well close to the top.</p> <p>After washing, blot the microwells on absorbent tissue. Follow wash protocol meticulously</p> <p>Use only Dengue NS1 Ag Microlisa wash solution.</p> <p>Incubate the plate in dark after addition of substrate.</p>
4. Poor reproducibility	<p>a) Washing problems.</p> <p>b) Uncalibrated pipettes or tips not well fitted, improper pipetting/ dispensing.</p> <p>c) Reagent & sera not at room temperature or not well mixed before use.</p> <p>d) Too long time for addition of controls, samples or reagents, Inconsistency in time intervals.</p> <p>e) Interference in optical pathway due to Air bubbles.</p>	<p>Check all 8 ports/ manifold for uniform flow of wash buffer. If there are blockage, clen the ports.</p> <p>Use only calibrated pipettes with well fitted tips & pipette carefully without bubbling.</p> <p>Equilibrate reagents to room temperature and mix thoroughly before use</p> <p>Develop consistent and uniform technique.</p> <p>Clean or dry the bottom of microwells, check for bubbles and repeat the readings.</p>

PROBLEM	POSSIBLE CAUSE	SOLUTION
5. False Positive	<p>a) Beside 3a, b, c, d, e & f incorrect interpretation and calculation of final results</p> <p>b) High incubator temperature, incorrect timing or pipetting</p> <p>c) Use of turbid/ lipaemic or hemolyzed sample.</p>	<p>Check the calculation part given in the insert and correctly interpret.</p> <p>Check incubator temperature, procedure & repeat assay.</p> <p>Centrifuge the sample at 5000 rpm for 30 minutes and re-run the test with clear sample.</p>
6. False Negative/ low O.D. of calibrator, Positive control & positive sample	<p>a) Inadequate addition of substrate/enzyme conjugate conc. solution.</p> <p>b) Kit expired, reagent of different kit used.</p> <p>c) White particles in working substrate solution.</p> <p>d) Uncalibrated pipettes, improper pipetting.</p> <p>e) Deterioration of calibrator and/or enzyme conjugate conc.</p> <p>f) Stop solution is added before 30 minutes. Reaction terminated before 30 minutes.</p> <p>g) O.D. taken at incorrect wavelength.</p> <p>h) Incorrect incubator temperature, timing or pipetting</p> <p>i) Kit deterioration</p> <p>j) Sample deterioration due to improper storage and / or microbial contamination.</p>	<p>Follow the procedure meticulously & repeat assay.</p> <p>Check the expiry of the kit before use.</p> <p>Discard the substrate and prepare the working substrate again in fresh tube.</p> <p>Use only calibrated pipettes with well fitted tips & pipette carefully without bubbling.</p> <p>Check storage of calibrator and/or Enzyme conjugate conc. They shall be stored at 2-8°C.</p> <p>Follow the test procedure meticulously.</p> <p>Read O.D. values at 450 nm and 630 nm.</p> <p>Check incubator temperature, procedure & repeat assay</p> <p>Check storage of kit and should be stored at 2-8°C.</p> <p>Store the sample at 2-8°C / -20°C as recommended in the specimen collection & handling.</p>

in vitro diagnostic Reagent, not for medicinal use

J. Mitra & Co. Pvt. Ltd.
A 180-181, Okhla Ind. Area, Ph-1, New Delhi-110 020, INDIA
Ph: +91-11-47130300, 47130500
e-mail: jmitra@jmitra.co.in Internet: www.jmitra.co.in