

DENGUE

A Continuing Global Threat.....



**Dengue can be
Diagnosed &
CONFIRMED by:**

ELISA RANGE:

NS1 ANTIGEN ELISA Test
Dengue NS1 Ag Microlisa

MAC ELISA Test
Dengue IgM Microlisa

GAC ELISA Test
Dengue IgG Microlisa

Tests also Available:

RAPID RANGE:

NS1 Antigen and IgM & IgG antibody Rapid Test
Dengue DAY1 Test

NS1 Antigen Rapid Test
Advantage Dengue NS1 Ag Card

POCT TEST:


NS1 Antigen (FINGER PRICK) Test
Dengue NS1 Antigen FP

FLUORESCENCE BASED TEST:

Dengue NS1 Antigen (FIA Test)
Dengue NS1 Ag Quanti Card

Dengue Antibody (FIA Test)
Dengue IgM Quanti Card
Dengue IgG Quanti Card

Dengue course of illness & approximate time-line of Primary and Secondary antibodies formation



Dengue is the most rapidly spreading mosquito-borne disease in the world. The key to handle dengue infection is its early recognition and understanding of the clinical problems during the different phases of the disease, leading to a rational approach to case management and good clinical outcome (Figure 2.1).

Dengue infection is the systematic and dynamic disease. After the incubation period, the illness begins abruptly and is followed by the 3 phases - febrile, critical and recovery (Figure 2.1).

NS1 ANTIGEN (Figure 4.1)

As can be observed from Figure 4.1, during early stage of disease, Antigen Detection can be used to diagnose the infection.

Detection of NS1 antigen is important for early and accurate diagnosis of dengue. NS1 Antigen can be detected from approximately Day 1 to Day 7 of fever setting in.

IgM ANTIBODIES (Figure 4.1)

Antibody response to infection differs according to the immune status of the host. In **primary infection**, IgM antibodies become detectable about 5-6 days after onset of disease. When the viremia declines, IgM level rises quickly to reach peak in about 2 weeks.

In **Secondary infection** IgM antibodies become detectable about 4-5 days and their levels are comparatively low.

IgG ANTIBODIES (Figure 4.1)

In **primary infection**, Antibody production of IgG will be at a lower level when compared to IgM. IgG antibodies are generally detectable at low levels in about 11-12 days and increase slowly and remain detectable after several months and probably even for life.

In **Secondary infection**, the IgG antibody level rises quickly reaching to peak in about 2 weeks after the onset of symptoms and may persists for years.

*source: Dengue Guidelines for Diagnosis, Treatment, Prevention & Control, New edition : 2009, WHO, Page No.: 25

**source: Dengue Guidelines for Diagnosis, Treatment, Prevention & Control, New edition : 2009, WHO, Page No.: 92

http://whqlibdoc.who.int/publications/2009/9789241547871_eng.pdf

NS1 Ag detection is the method of choice being the direct method and have a balance relationship, between the ease of use or accessibility of a diagnostic method and the confidence in the results of the test (Figure 4.2).


*** source: Dengue Guidelines for Diagnosis, Treatment, Prevention & Control, New edition: 2009, WHO, Page No.: 92

http://whqlibdoc.who.int/publications/2009/9789241547871_eng.pdf

For Educational Purpose only.




FIGURE 4.2 COMPARISON OF DIAGNOSTIC TESTS ACCORDING TO THEIR ACCESSIBILITY AND CONFIDENCE ***



Dengue NS1 Ag Microlisa

Microwell ELISA Test for the Detection of Dengue NS1 Antigen in Human Serum/ Plasma

TEST PRINCIPLE



- Direct test method for Detecting infection with Dengue NS1 Ag
- Excellent sensitivity against all the 4 serotypes of Dengue Virus Antigen
- Significant reduction of Window period
- Based on Direct Sandwich ELISA
- No sample dilution is required
- Long Shelf Life: 24 months at 2-8°C
- Color Coded reagent
- Convenient Packsize: 96 Tests

TEST PROCEDURE

Add 50 µl diluent in each well. Add 50 µl NC in A1 well, 50µl calibrator in B1, C1 & D1well, add 50 µl PC in E1 well, and add 50 µl sample to F1 well onwards																																																					
Prepare working conjugate solution	<table border="1"> <thead> <tr> <th>No. of strip</th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> <th>6</th> <th>7</th> <th>8</th> <th>9</th> <th>10</th> <th>11</th> <th>12</th> </tr> </thead> <tbody> <tr> <td>Enzyme</td> <td>20</td> <td>40</td> <td>60</td> <td>80</td> <td>100</td> <td>120</td> <td>140</td> <td>160</td> <td>180</td> <td>200</td> <td>220</td> <td>240</td> </tr> <tr> <td>Conjugate (µl)</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Diluent (ml)</td> <td>1</td> <td>2</td> <td>3</td> <td>4</td> <td>5</td> <td>6</td> <td>7</td> <td>8</td> <td>9</td> <td>10</td> <td>11</td> <td>12</td> </tr> </tbody> </table>	No. of strip	1	2	3	4	5	6	7	8	9	10	11	12	Enzyme	20	40	60	80	100	120	140	160	180	200	220	240	Conjugate (µl)													Diluent (ml)	1	2	3	4	5	6	7	8	9	10	11	12
No. of strip	1	2	3	4	5	6	7	8	9	10	11	12																																									
Enzyme	20	40	60	80	100	120	140	160	180	200	220	240																																									
Conjugate (µl)																																																					
Diluent (ml)	1	2	3	4	5	6	7	8	9	10	11	12																																									
Add 100 µl working conjugate in each well																																																					
Cover the plate and incubate for 90 min. at 37°C																																																					
Wash 6 cycles																																																					
Prepare working substrate solution	<table border="1"> <thead> <tr> <th>No. of strip</th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> <th>6</th> <th>7</th> <th>8</th> <th>9</th> <th>10</th> <th>11</th> <th>12</th> </tr> </thead> <tbody> <tr> <td>TMB Substrate (ml)</td> <td>0.8</td> <td>1.6</td> <td>2.4</td> <td>3.2</td> <td>4.0</td> <td>4.8</td> <td>5.6</td> <td>6.4</td> <td>7.2</td> <td>8.0</td> <td>8.8</td> <td>9.6</td> </tr> <tr> <td>TMB Diluent (ml)</td> <td>0.8</td> <td>1.6</td> <td>2.4</td> <td>3.2</td> <td>4.0</td> <td>4.8</td> <td>5.6</td> <td>6.4</td> <td>7.2</td> <td>8.0</td> <td>8.8</td> <td>9.6</td> </tr> </tbody> </table>	No. of strip	1	2	3	4	5	6	7	8	9	10	11	12	TMB Substrate (ml)	0.8	1.6	2.4	3.2	4.0	4.8	5.6	6.4	7.2	8.0	8.8	9.6	TMB Diluent (ml)	0.8	1.6	2.4	3.2	4.0	4.8	5.6	6.4	7.2	8.0	8.8	9.6													
No. of strip	1	2	3	4	5	6	7	8	9	10	11	12																																									
TMB Substrate (ml)	0.8	1.6	2.4	3.2	4.0	4.8	5.6	6.4	7.2	8.0	8.8	9.6																																									
TMB Diluent (ml)	0.8	1.6	2.4	3.2	4.0	4.8	5.6	6.4	7.2	8.0	8.8	9.6																																									
Add 150 µl working substrate in each well																																																					
Incubate in dark for 30 min at room temperature. (20°C - 30°C)																																																					
Add 100µl Stop Solution and read result at 450 nm/630 nm.																																																					



Sensitivity: 99.5%
Specificity: 100%

Dengue IgM Microlisa











ELISA TEST

Dengue IgG Microlisa

MAC ELISA (Capture Principle) Test for the Detection of Dengue IgM Antibodies in Human Serum/ Plasma


GAC ELISA (Capture Principle) Test for the Detection of Dengue IgG Antibodies in Human Serum/ Plasma




TEST PROCEDURE

Dilute sample (1:100) sample 10 µl + sample diluent 1000µl	
Add 100 µl control, calibrator & diluted sample in respective wells.	
Cover the plate and incubate for 60 min. at 37°C	
Wash (5 cycles)	
Add 100 µl working conjugate in each well	
Cover the plate and incubate for 60 min. at 37°C	
Wash (5 cycles)	
Prepare working substrate solution	No. of strip 1 2 3 4 5 6 7 8 9 10 11 12 TMB 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0 5.5 6.0 Substrate (ml). TMB 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0 5.5 6.0 Diluent (ml).
Add 100 µl working substrate in each well	
Incubate in dark for 30 min at room temperature.	
Add 50 µl Stop Solution and read result at 450 nm/630 nm.	

* The Preparation of working substrate solution mentioned above in the table is suggestive only, however working substrate solution can be prepared in any amount proportionately as per requirement of tests to be performed (from 5 tests to 96 tests as 3 calibrator, 1 negative control and 1 positive control are to be included in each run).

TEST PRINCIPLE



-  Anti-Human IgM/IgG antibody
-  IgG antibody in the sample
-  Substrate TMB
-  IgM antibody in the sample
-  Dengue Antigen conjugated with HRPO
-  Stop Solution



Sensitivity: 99.13%

Specificity: 99.84%

- Principle : "Advance Capture" ELISA
- Enzyme Conjugate (Ready to use), liquid stable for 18 months, no wastage
- Excellent Sensitivity against all 4 serotypes of Dengue
- High signal to noise ratio
- Color Coded reagents to monitor procedural steps
- Long shelf Life: 18 months at 2-8°C
- Unique Break-away wells
- Convenient Packsize: 96 Tests



Sensitivity: 98.66%

Specificity: 99.93%

RAPID TEST



Dengue NS1 Ag:

Sensitivity: 96% Specificity: 98%

Dengue IgM/IgG antibody:

* Patented Device for Dengue NS1 Ag
** Patented Device for Dengue IgM & IgG Antibody Card

- First-Line Testing Kit for Dengue diagnosis
- Rapid visual test for Dengue NS1 Antigen and IgM & IgG Antibodies detection
- Diagnosis of both Primary & Secondary Infection
- Detects all 4 serotypes of Dengue Virus
- Highly Sensitive & Highly Specific
- Long Shelf Life 30 months at 2-30°C
- Convenient Packsize: 10 Tests, 30 Tests & 50 Tests

“Why Dengue DAY1 Test, Antigen & Antibody (Combi test) is more Reliable in Dengue Diagnosis”

Study highlights that **using dengue NS1 antigen detection in combination with IgM and IgG serology can significantly increase the sensitivity of acute dengue diagnosis** and extends the possible window of detection to include very early acute samples and enhances the clinical utility of rapid immunochromatographic testing for dengue.

Source : RESEARCH ARTICLE: The Diagnostic Sensitivity of Dengue Rapid Test Assays Is Significantly Enhanced by Using a Combined Antigen and Antibody Testing Approach
www.plosntds.org/article/info%3Adoi%2F10.1371%2Fjournal.pntd.0001199

POCT (Point of Care) TEST

**Ideal test for
Physician's Clinic
&
Nursing Homes
for immediate
results.**



DENGUE NS1 ANTIGEN FP

Finger Prick quick visual test for the detection of Dengue NS1 Antigen in Human Whole Blood/ Serum/ Plasma

- **Point of Care Test (POCT)**
- Finger Prick quick visual test for detection of Dengue NS1 Antigen
- Results within 20 minutes
- From DAY1 of the fever, It detects all 4 serotypes of Dengue Virus
- Highly Sensitive & Highly Specific
- Individually pouched test device
- Longer Shelf Life: 24 Months At 2-30°C
- Patented Product: Patent No.: **262648**
- Prompt Service Backup
- Convenient Packsize: 10 Tests

**Advanced Technology
Accurate Results
Affordable Price**

Developed in partnership
with HTIC,
IIT Madras, India



**DENGUE NS1 Ag
Quanti Card**

Fluorescence immunoassay for qualitative
detection of Dengue NS1 Antigen
in Human Serum/ Plasma

- ▶ **Principle: Based on most advanced fluorescence immunoassay method**
- ▶ **Individually pouched test cartridge, each cartridge is precision engineered with embeded lot data.**
- ▶ **Results within 30 Minutes**
- ▶ **Sensitivity: 100%* & Specificity: 98.38%* (In comparison with ELISA Test)**
- ▶ **Long Shelf Life: 24 Months at 2-8°C**
- ▶ **Available Packsize: 24 Tests**

* Source: <https://doi.org/10.20546/ijcmas.2017.611.061>, Table-3, Page- 510



**DENGUE IgM
Quanti Card**

Fluorescence immunoassay for the qualitative
measurement of Dengue IgM Antibodies
in Human Serum/ Plasma

- ▶ **Principle: Based on most advanced fluorescence immunoassay method**
- ▶ **Individually pouched test cartridge, each cartridge is precision engineered with embeded lot data.**
- ▶ **Results within 40 Minutes**
- ▶ **Sensitivity: 96.72% & Specificity: 97.45% (In comparison with ELISA Test: Dengue IgM Microlisa)**
- ▶ **Long Shelf Life: 24 Months at 2-8°C**
- ▶ **Available Packsize: 24 Tests**



**DENGUE IgG
Quanti Card**

Fluorescence immunoassay for the qualitative
measurement of Dengue IgG Antibodies
in Human Serum/ Plasma

- ▶ **Principle: Based on most advanced fluorescence immunoassay method**
- ▶ **Individually pouched test cartridge, each cartridge is precision engineered with embeded lot data.**
- ▶ **Results within 40 Minutes**
- ▶ **Sensitivity: 97.23% & Specificity: 97.36% (In comparison with ELISA Test: Dengue IgG Microlisa)**
- ▶ **Long Shelf Life: 24 Months at 2-8°C**
- ▶ **Available Packsize: 24 Tests**

Also Connect with us on:



ISO 13485:2016
ICMED 13485



For further enquiries, Please contact:
J. Mitra & Co. Pvt. Ltd.

A-180-181, Okhla Indl. Area, Phase-1, New Delhi - 110 020, INDIA

Tel: +91-11-47130300, 47130500

E-mail: jmitra@jmitra.co.in Website : www.jmitra.co.in