


MALARIA

re-emerging as world No.1 infectious Killer

Malaria is the most common infectious tropical disease and an enormous public health problem. It is caused by the parasite called plasmodium and is passed to human by mosquito (female anopheles). It is a disease that can be treated in just 2 days, yet it can cause fatal complication if the *diagnosis is delayed*.


Antigen Test Made Malarial Detection Just Early, Easy & Accurate.....

- Advantage Mal Card
- Advantage Malaria Pan+Pf Card
- Advantage Pf Malaria Card
- Advantage Malaria Card
- Advantage Pan Malaria Card



Combining Reliability with Convenience through Rapid Antigen Test...

Lets know about J. Mitra Malaria principle.....



PLASMODIUM SPECIES DETECTION BY DIFFERENT TEST

Product Name	Test Line	PAN	P.v	Pf	
				HRP-2	pLDH
Advantage Mal Card		✓	—	—	✓
Advantage Malaria Pan+Pf Card		✓	—	✓	—
Advantage Pf Malaria Card		—	—	✓	—
Advantage Malaria Card		—	✓	✓	—
Advantage Pan Malaria Card		✓	—	—	—

Advantage MAL CARD

For detection of *P. falciparum* - pLDH & Plasmodium Species (*P. falciparum*/*P. vivax*/*P. malariae*/*P. ovale* - pLDH) in Human Whole Blood



Sensitivity: **100%** Specificity: **99.95%**
In-house evaluation

- Individually pouched Single Tests - Practical & Safe
- Analytical Sensitivity of 100-200 parasite per μ l blood
- See through Device
- Longer Shelf Life 30 months at 4-30°C
- Easy to Interpret Colour Bands
- Result within 20 minutes
- Convenient Pack size: 50 Tests & 100 Tests
- Panel Detection Score: **97.1%***(P.v) Specificity: **99.60%***
*As per WHO, Geneva evaluation

Advantage Malaria PAN + P.f CARD

For detection of *P. falciparum* - HRP-2 & Plasmodium Species (*P. falciparum*/*P. vivax*/*P. malariae*/*P. ovale* - pLDH) in Human Whole Blood



- Individually pouched Single Tests - Practical & Safe
- Analytical Sensitivity of 50-100 parasite per μ l blood
- See through Device
- Longer Shelf Life 30 months at 4-30°C
- Easy to Interpret Colour Bands
- Result within 20 minutes
- Convenient Pack size: 50 Tests & 100 Tests
- Panel Detection Score: **100%*** (P.v & Pf) Specificity: **100%***
*As per WHO, Geneva evaluation

Advantage P.f MALARIA CARD

For detection of *P. falciparum* - HRP-2 in Human Whole Blood



Sensitivity: **100%** Specificity: **99.80%**
In-house evaluation

- Individually pouched Single Tests - Practical & Safe
- Analytical Sensitivity of ≥ 50 parasite per μ l blood
- See through Device
- Longer Shelf Life 30 months at 4-30°C
- Easy to Interpret Colour Bands
- Result within 20 minutes
- Convenient Pack size: 25 Tests
- Panel Detection Score: **99%*** Specificity: **100%***
*As per WHO, Geneva evaluation

Advantage MALARIA CARD

For detection of *P. falciparum* - HRP-2 & Plasmodium Species (*P. vivax* - pLDH) in Human Whole Blood



Sensitivity: **100%** Specificity: **99.77%**
In-house evaluation

- Individually pouched Single Tests - Practical & Safe
- Analytical Sensitivity of 50-200 parasite per μ l blood
- See through Device
- Longer Shelf Life 30 months at 4-30°C
- Easy to Interpret Colour Bands
- Result within 20 minutes
- Convenient Pack size: 25 Tests & 50 Tests

Advantage PAN MALARIA CARD


For detection of Plasmodium Species (*P. falciparum*/*P. vivax*/*P. malariae*/*P. ovale* - pLDH) in Human Whole Blood



Sensitivity: **100%** Specificity: **99.8%**
In-house evaluation

- Individually pouched Single Tests - Practical & Safe
- Analytical Sensitivity of 100 parasite per μ l blood
- See through Device
- Longer Shelf Life 30 months at 4-30°C
- Easy to Interpret Colour Bands
- Result within 20 minutes
- Convenient Pack size: 50 Tests & 100 Tests
- Panel Detection Score: **100%*** (P.v) Specificity: **99.6%***
*As per WHO, Geneva evaluation

TEST PROCEDURE



Add 4 μ l of the anti-coagulated blood sample using sample dropper or micropipette (or finger prick sample) to the sample well 'S'. Make sure that the blood from the sample dropper has been completely transferred to the sample pad.

Add 3 drops of Assay Buffer to buffer well 'B' and allow the reaction to occur for 20 minutes.

Read Result in 20 minutes.

Important Note: Discard the DEVICE immediately after reading results as it is potentially infectious.

*Products Performance as per WHO, Geneva and FIND, Geneva

in 2013, WHO, Geneva and FIND (Foundation for Innovative New Diagnostics), Geneva launched an evaluation program to assess the performance of commercially available malaria RDTs and their results are:

Figure 13: Phase-2 *P. vivax* panel detection score and positivity rate at 200 parasites/ μ L





Figure 15: Phase-2 Plasmodium spp. (pan or *P. vivax*/Pvom test line) false-positive rate against clean-negatives



Malaria Rapid Diagnostic Test Performance, Results of WHO product testing of malaria RDTs: Round 5 (2013), Page No.: 46
Website: http://apps.who.int/iris/bitstream/10665/128678/1/9789241507554_eng.pdf

* Note : The above information is provided for the scientific community, It is not for commercial or promotional purpose.

Also Connect with us on:



For further enquiries, Please contact:

J. Mitra & Co. Pvt. Ltd.

A 180-181, Okhla Indl. Area, Phase-1, New Delhi - 110 020, INDIA

Tel: +91-11-47130300, 47130500

E-mail: jmitra@jmitra.co.in Website: www.jmitra.co.in