



- Keep the D-DIMER conjugate vial undisturbed for 1 minute & then vortex it well for 30 seconds.
- Immediately load 50µl of above sample reagent mixture using micropipette to the sample well of the Cartridge. **Care should be taken to avoid any spillage on the QR-Code and test result window.**
- Leave the sample loaded cartridge at room temperature for 20 minutes.
- After 20 minutes, gently add 2 drops of assay buffer in the buffer well of the test cartridge and leave the cartridge for another 10 minutes at room temperature.
- In the meantime, enter the patient's details in the iQuant analyzer testing window and select the D-DIMER test from the pop down menu in the testing window of the iQuant analyzer.
- After the incubation period is over, Insert the test cartridge into the iQuant Analyzer with arrow (←) marked side on the top of cartridge facing towards the analyzer and press RUN button. Note down the value displayed on the screen of iQuant Analyzer.

**Discard the D-Dimer Quanti Card immediately after reading results considering it to be potentially infectious.**

<b>D-DIMER REFERENCE RANGE</b>	: <500 ng/ml.
<b>MEASURING RANGE</b>	: 50-10,000 ng/ml.
<b>DETECTION LIMIT</b>	: 50 ng/ml.

### PERFORMANCE CHARACTERISTICS OF D-DIMER QUANTI CARD

#### 1. Precision

**Intra-Assay:** Within-run and between-run precision have been determined by testing 10 replicates of 3 different samples with D-Dimer concentration: 238.4ng/ml, 3224.5ng/ml and 7512.6ng/ml on the three different lots. The C.V (%) for all 3 samples is ≤ 10%.

**Inter-Assay:** The inter-assays were performed with 10 replicates of 3 different samples with D-Dimer concentration: 238.4ng/ml, 3224.5ng/ml and 7512.6ng/ml on three different lots on 10 sequential days. The C.V (%) for all 3 samples is ≤ 10%.

#### 2. Accuracy

The accuracy of D-Dimer Quanti Card was checked with 75 clinical specimens, D-Dimer concentrations of 75 samples were compared with commercially available kit. All the samples are run in duplicate on D-Dimer Quanti Card. The following results were obtained:

Slope	: 0.9842
Y-Intercept	: 80.076
R <sup>2</sup>	: 0.996

#### 3. Specificity

There was no significant interference with the D-Dimer measurement when other biomolecules such as Bilirubin (20mg/dL), Triglyceride (1000mg/dL), Glucose (500mg/dL), Haemoglobin (500mg/dL) and Albumin (10g/dL) were added to the test specimen with much higher level than in normal blood.

### LIMITATIONS AND INTERFERENCES

- The test procedure, precautions and interpretation of results for this test must be strictly followed.
- As with all diagnostic tests, the test result must always be correlated with clinical finding and laboratory data available.
- Any modification to the above procedure and / or use of other reagents will invalidate the test procedure.
- The presence of additional substances in blood samples may interfere with product performance and may cause erroneous results.

### LIMITED EXPRESSED WARRANTY DISCLAIMER

The manufacturer limits the warranty to the test kit, as much as that the test kit will function as an in vitro diagnostic assay within the limitations and specifications as described in the product instruction for use, when used strictly in accordance with the instructions contained therein. The manufacturer disclaims any warranty expressed or implied including such expressed or implied warranty with respect to merchantability, fitness for use or implied utility for any purpose. The manufacturer's liability is limited to either replacement of the product or refund of the purchase price of the product and in no case liable to claim of any kind for an amount greater than the purchase price of the goods in respect of which damages are likely to be claimed. The manufacturer shall not be liable to the purchaser or third parties for any injury, damage or economic loss, howsoever caused by the product in the use or in the application there of.

### TROUBLE SHOOTING

PROBLEM	POSSIBLE CAUSE	SOLUTION
1. Unexpected low Test Result	a) Hemolyzed sample b) Use of serum sample for testing c) Insufficient volume of sample used.	Repeat the test using fresh Plasma samples.

d) Fresh Plasma sample not used.

e) Reading has been taken at less than the prescribed time.

f) Reagents used were too cold and were not brought to Room Temperature RT.

g) D-Dimer Conjugate is continuously exposed to light for more than 8 hours.

h) Expired Test kit used.

i) Improper i.e. less volume of reaction mixture applied to sample well of cartridge.

a) High amount of plasma sample used.

b) Results read beyond 10 minutes after adding assay buffer

c) Improper i.e. high volume of reaction mixture applied to sample well of cartridge

Read the results at prescribed time only.

Bring the whole test kit RT before testing

Always reseal the unused D-Dimer Conjugate Vials with clamp & rod after use. Store at 2-8 °C.

Repeat the test using a new test kit that has not passed the expiration date.

Use appropriate volume of reaction mixture using calibrated pipette.

Retest using 100µl sample volume

Read the test result at 10 minutes.

Use 50µl of reaction mixture using calibrated pipette.

2. Unexpected high test result

### REFERENCES

- Yao Y, Cao J, Wang Q, et. al., D-dimer as a biomarker for disease severity and mortality in COVID-19 patients: a case control study. Journal of Intensive Care 2020; 8, 49-60.
- Favressea J, Lippig G, Roy PM, et. al., D-dimer: Preanalytical, analytical, postanalytical variables, and clinical applications. Critical Reviews in Clinical Laboratory Sciences 2018; 55(8) 548-577.
- Adam SS, Key NS, Greenberg CS. D-dimer antigen: current concepts and future prospects. Blood 2009; 113, 2878-2887.

in vitro diagnostic reagent, not for medicinal use

**J. Mitra & Co. Pvt. Ltd.**  
 A 180-181, Okhla Ind. Area, Ph-1, New Delhi-110 020, INDIA  
 Ph: +91-11-47130300, 47130500  
 e-mail: jmitra@jmitra.co.in Internet: www.jmitra.co.in