


- Add whole blood sample using the sample dropper/ micropipette onto the sample pad in the sample well 'S'. fig. (g)

NOTE: Press the tip of the dropper onto the sample pad in the sample well "S" to ensure that the complete volume of whole blood has been transferred to the strip.


- Add 3 drops of the Assay Buffer in the buffer well 'B'. fig. (h). **Ensure FREE FALLING OF DROPS on the membrane, holding the vial/dropper vertically for proper volume. Screw cap the vial after use.**




- Allow the reaction to occur for 20 minutes.
- Read the results at 20 minutes. **Do not read the result after 30 minutes. Reading beyond prescribed time may give false results.**
- Discard the Advantage Malaria Pan + Pf Card immediately after reading results at 20 minutes as it is potentially infectious.

INTERPRETATION OF THE RESULTS


REACTIVE



As shown in fig. (i), appearance of three purplish pink coloured lines one each in Pf. region (F), Pan region (P) & Control region (C) indicates that the sample is reactive for P. falciparum or mixed infection of Pf and Pv (or P. malariae, Povale).




As shown in fig. (j) appearance of two purplish pink coloured line one each at P & C region only indicates that the sample is reactive for P. vivax / P. malariae / P. ovale only.



As shown in fig. (k) appearance of two purplish pink coloured line one each at F & C region only indicates that the sample is reactive for P. falciparum only. A difference of intensity in colour may occur between both the test lines ('P' & 'F') and between the test lines & control line depending on the concentration of pLDH & HRP-2 in the sample but this does not affect the interpretation of the results.



Depending on the concentration of pLDH & HRP-2 positive results may be observed within 60 seconds. However, to confirm a negative result the test result should be read only at 20 minutes. **Consider a faint test line also as a positive result.**

NON-REACTIVE



As shown in fig. (l), appearance of only one purplish pink coloured line at Control (C) region indicates that the sample is non-reactive for P. falciparum and other Plasmodium Species (P. vivax / P. malariae / P. ovale)

INVALID



The test is invalid, if no control line appears after the completion of test, either with clear background or with complete pinkish/ purplish background fig. (m, n, o & p). **The test should be repeated using a new card.**

LIMITATIONS AND INTERFERENCES

- The test procedure, precautions and interpretation of results for this test must be strictly followed.
- This is only a screening test. All reactive samples must be confirmed with microscopy. As with all diagnostic tests, the test result must always be correlated with clinical finding. The results should be reported only after complying the mentioned procedure.
- Though the test is accurate in detecting HRP-2 specific to P. falciparum or pLDH specific to Plasmodium species (P. falciparum/ P. vivax/ P. malariae/ P. Ovale), a low incidence of false results can occur. All reactive test results are to be interpreted within the epidemiological, clinical and therapeutic context. When it seems indicated, the parasitological techniques of reference should be considered (microscopic examination of the thick smear and thin blood films).
- Any modification to the test procedure and/ or use of reagents other than provided with the test kit will lead to invalid and/ or false test result.
- Since the HRP-2 persists for upto a fortnight even after successful anti-malarial treatment, a positive test result does not indicate a failed therapeutic response.
- In P. falciparum malaria infection, HRP-2 is not secreted in gametogony stage. Hence, in 'carriers', the HRP-2 test line (F) may be absent.
- The possibility of resistant strain of malaria should always be considered if the reaction of the test remains positive with the same intensity after 5-10 days post treatment.
- Patient with rheumatoid factors, anti-nuclear antibody or dengue may give false positive results.
- The product performance will be hampered or degraded if test kit is stored at low temperature (<2°C) or high temperature (>30°C).
- False negative results may be obtained in following conditions:

i) Parasite density/ antigen concentration is below the detection limit of the test or analyte detected are not present during the stage of disease in which specimen is collected.

ii) No production of HRP-2 antigen in the specimen due to deletion of HRP-2 gene.

- In case of very faint or doubt for test band (F and/ or P), the test should be repeated using fresh device.
- Repeat the test in case of strong clinical evidence of malaria using fresh device.

PERFORMANCE CHARACTERISTICS OF ADVANTAGE MALARIA PAN + Pf CARD

(i) WHO Evaluation:

The ADVANTAGE MALARIA PAN + Pf CARD test kit has been evaluated by WHO, Geneva using a panel of wild & cultured malaria positive sample and the results obtained are as follows :

	Panel Detection Score*		Specificity*
	Pf	Pv	
200 Parasites	84%	100%	—
2000 Parasites	100%	100%	—
Negative Samples	—	—	100%

*Reference: Malaria Rapid Diagnostic Test Performance: Results of WHO product testing of malaria RDTs: Round 5 (2013), Page:38.

Note : The above information is provided for the scientific community, It is not for commercial or promotional purpose.

(ii) In-house Evaluation:

Analytical Sensitivity: The test can detect parasitemia levels of ≥ 50 parasites per μ l of blood for P. falciparum (HRP-2), ≥ 100 parasites per μ l of blood for P. falciparum (pLDH) & P. vivax.

The ADVANTAGE MALARIA PAN + Pf CARD has been evaluated in-house with malaria positive and negative clinical whole blood samples and compared with microscopic examination. The evaluation also included cross-reacting samples; Dengue, Rheumatoid factor, Leptospira, HIV, HCV, HBV, M. tuberculosis, Syphilis, Brucella, Scrub typhus positive samples. The results obtained are as follows:

Sample	Total no. of samples tested	ADV. MALARIA PAN + Pf CARD		Sensitivity (%)	Specificity (%)
		Positive	Negative		
Malaria Negative	2100	1	2099	-	99.95
Cross-reacting sample	64	0	64	-	100
P. falciparum Positive	58	58	0	100	-
P. vivax Positive	105	105	0	-	100

Precision: Within-run and between-run precisions have been determined by testing 10 replicates of six specimens : two negative, two weak positive and two moderate positive. The C.V (%) of negative, weak positive and moderate positive samples were within 10% of the time.

LIMITED EXPRESSED WARRANTY DISCLAIMER

The manufacturer limits the warranty to the test kit, as much as that the test kit will function as an in vitro diagnostic assay within the limitations and specifications as described in the product instruction-manual, when used strictly in accordance with the instructions contained therein. The manufacturer disclaims any warranty expressed or implied including such expressed or implied warranty with respect to merchantability, fitness for use or implied utility for any purpose. The manufacturer's liability is limited to either replacement of the product or refund of the purchase price of the product and in no case liable to claim of any kind for an amount greater than the purchase price of the goods in respect of which damages are likely to be claimed. The manufacturer shall not be liable to the purchaser or third parties for any injury, damage or economic loss, howsoever caused by the product in the use or in the application thereof.

BIBLIOGRAPHY OF SUGGESTED READING

- Jamshaid Iqbal; etal. Journal of Clinical Microbiology, 37(1), 1999, 3644-3646.
- World Health Organization 2000. New perspectives malaria diagnosis. World Health Organization, Geneva, Switzerland.
- Piper R, etal. Am J. Trop. Med. Hyg., 60(1), 1999, 109-118.
- Moody A, Clinical Microbiological Review, 15(1), 2002, 66-78.
- Malaria Rapid Diagnostic Test Performance: Results of WHO product testing of malaria RDTs: Round 5 (2013)

in vitro diagnostic reagent, not for medicinal use

J. MITRA & CO. PVT. LTD.

A 180-181, Okhla Indl. Area, Phase-1, New Delhi-110 020, INDIA

Ph.: +91-11-47130300, 47130500

e-mail: jmitra@jmitra.co.in Internet: www.jmitra.co.in