

DENGUE NS1 Ag MICROLISA

Microwell ELISA Test for the Detection of Dengue NS1 Antigen in Human Serum/Plasma

1. INTRODUCTION

Dengue virus is a flavivirus found largely in areas of the tropic and sub-tropics. There are four distinct but antigenically related serotypes of dengue viruses, and transmission is by mosquito, principally Aedes aegypti and Aedes albopictus.

The mosquito-borne dengue viruses (serotype 1-4) cause dengue fever, a severe flu-like illness. The disease is prevalent in third world tropical regions and spreading to sub-tropical developed countries - including the United States. WHO estimates that 50-80 million cases of dengue fever occur worldwide each year, including a potentially deadly form of the disease called dengue haemorrhagic fever (DHF) and dengue shock syndrome (DSS). Primary infection with dengue virus results in a self-limiting disease characterized by mild to high fever lasting 3 to 7 days, severe headache with pain behind the eyes, muscle and joint pain, rash and vomiting. Secondary infection is the more common form of the disease in many parts of Southeast Asia and South America. This form of the disease is more serious and can result in DHF and DSS. The major clinical symptoms can include high fever, haemorrhagic evets, and circulatory failure, and the fatality rate can be as high as 40%. Early diagnosis of DSS is particularly important, as patients may die within 12 to 24 hours if appropriate treatment is not administered.

Primary dengue virus infection is characterized by elevations in specific NS1 antigen levels 0 to 9 days after the onset of symptoms; this generally persists up to 15 days. Consequently, testing for the dengue NS1 antigen via ELISA is a highly effective diagnostic tool because it enables early diagnosis prior to seroconversion. Earlier diagnosis of Dengue reduces risk of complication such as DHF or DSS, especially in countries where dengue is endemic. Unlike early antigens **IgM antibodies** do not become detectable until 3 to 5 days after symptom onset. Secondary dengue virus infections are characterized by rapidly rising **IgG levels**, which peak within a 6-to-15-day detection window and may be accompanied by elevated IgM. The J. Mitra Dengue IgG Microlisa is specifically engineered to detect these high, diagnostic thresholds of IgG antibodies. Therefore, to ensure maximum diagnostic accuracy across various clinical stages, a combined testing protocol using both the **J. Mitra Dengue IgM Microlisa** and the **Dengue IgG Microlisa** is highly recommended.

2. INTENDED USE

DENGUE NS1 Ag MICROLISA is an in vitro diagnostic medical device intended for the qualitative detection of Dengue Virus Non-Structural Protein 1 (NS1) antigen in human serum or plasma specimens. The assay is intended to aid in the early diagnosis of acute dengue virus infection caused by Dengue virus serotypes DEN-1, DEN-2, DEN-3, and DEN-4 in individuals suspected of dengue infection based on clinical symptoms and epidemiological risk factors.

The device is a manual microplate ELISA-based assay intended for professional laboratory use by trained healthcare professionals. The assay provides qualitative results and is intended to be used in conjunction with clinical findings, patient history, and other laboratory diagnostic information. The test is not intended as a standalone confirmatory assay.

3. PRINCIPLE

DENGUE NS1 Ag MICROLISA is a solid phase enzyme linked immunosorbent assay (ELISA) based on the "Direct Sandwich" principle. The microwells are coated with Anti-dengue NS1 antibodies with high reactivity for Dengue NS1 Ag. The samples are added in the wells followed by addition of enzyme conjugate (monoclonal anti-dengue NS1 antibodies linked to Horseradish Peroxidase (HRPO)). A sandwich complex is formed in the well wherein dengue NS1 (from serum sample) is "trapped" or "sandwiched" between the antibody and antibody HRPO conjugate. Unbound conjugate is then washed off with wash buffer. The amount of bound peroxidase is proportional to the concentration of dengue NS1 antigen present in the sample. Upon addition of the substrate buffer and chromogen, a blue colour develops. The intensity of developed blue colour is proportional to the concentration of dengue NS1 antigen in sample. To limit the enzyme-substrate reaction, stop solution is added and a yellow colour develops which is finally read at 450nm spectrophotometrically.

4. DESCRIPTION OF SYMBOLS USED

The following are graphical symbols used in or found on J. Mitra diagnostic products and packing. These symbols are the most common ones appearing on medical devices and their packing. They are explained in more detail in European Standard EN ISO 15223-1:2021.

	Manufactured By		In vitro diagnostic medical device
	No. of tests		Instruction for use
	Lot Number Batch Number		Temperature Limitation
	Manufacturing Date		Caution, see instruction for use
	Expiry Date		Catalogue Number
	Keep away from sunlight		Do not use if package is damaged
	Contains biological Material of Human Origin		Contains biological Material of Animal Origin
	Country of Manufacture		Keep Dry

5. KIT PRESENTATION

- 96 Tests

6. KIT & ITS COMPONENTS

Microwells	Microwells coated with anti-Dengue NS1 antibodies packed in a sealed pouch with desiccant.
Diluent	Buffer containing protein stabilizers & antimicrobial agents as preservative and to be used for Sample & Conjugate dilution.
Enzyme Conjugate Concentrate (50X)	Containing Monoclonal Anti-Dengue NS1 linked to horseradish peroxidase with protein stabilizers.
Wash Buffer Concentrate (25x)	Concentrated Phosphate buffer with surfactant.
TMB Substrate	To be diluted with TMB diluent before use.
TMB Diluent	Buffer solution containing H ₂ O ₂ with preservative
Control	Normal human serum negative for Dengue NS1 antigen with preservative.
Control	Recombinant Dengue NS1 antigen, with preservative.
Calibrator	Recombinant Dengue NS1 antigen, with preservative.
Stop Solution	Ready to use, 1N H ₂ SO ₄ .
Plate Sealers	Adhesive sheet to cover the microwell during incubation.

PROBLEM	POSSIBLE CAUSE	SOLUTION
	i) Kit deterioration	Check storage of kit and it should be stored at 2-8°C.
	j) Sample deterioration due to improper storage and / or microbial contamination.	Store the sample at 2-8°C / -20°C as recommended in the specimen collection & handling.



J. MITRA & CO. PVT. LTD.
A 180-181, Okhla Indl. Area, Phase-1, New Delhi-110 020, INDIA
Ph: + 91-11-47130300, 47130500
e-mail: jmitra@jmitra.co.in Internet: www.jmitra.co.in

in vitro diagnostic Reagent, not for medicinal use

VER-03
Rev. Date: Apr.-26

ANVEDA/007

R-00

PROBLEM	POSSIBLE CAUSE	SOLUTION
	f) Incorrect (low) incubator temperature, timing or pipetting g) Improper preparation of wash buffer, error of dilution, improper mixing of reagents. h) Kit deterioration	Check incubator temperature, procedure & repeat assay Check procedure & repeat assay Check storage of the kit and should be store at 2-8°C.
2. High O.D. value of Negative control	a) Plate not stopped after 20 minutes of adding stop solution b) Same microtip used for calibrator, positive and negative controls c) Nonspecific attachment/ binding of other reagent	Follow the procedure meticulously & repeat assay. Change micropipette tips while addition of negative/ positive control or calibrator If plates get scratches/ aberrations during washing, non specific proteins may bind while addition of next step.
3. Too much colour in all wells of the plate (high background)	a) Contaminated substrate b) Contaminated washing solution (1X). Poor quality of water used for diluting wash buffer conc. c) Over incubation of substrate and delay in addition of stop solution. d) Insufficient washing. i) Washing not consistent due to blockage of probes ii) Filling volume not sufficient. iii) Insufficient no. of wash cycles. iv) Contaminated wash device e) Use of wash solution from other manufacturer. f) Working substrate not protected from light	Check substrate (TMB Diluent) it should be colourless. If blue in colour then discard and use clean disposable container. Check the container and quality of water used for dilution. Use of distilled water is preferred. Follow the procedure meticulously. Check wash device and clean probes of manifold. fill the well close to the top. After washing, blot the microwells on absorbent tissue. Follow wash protocol meticulously Use only Dengue NS1 Ag Microlisa wash solution. Incubate the plate in dark after addition of substrate.
4. Poor reproducibility	a) Washing problems. b) Uncalibrated pipettes or tips not well fitted, improper pipetting/ dispensing. c) Reagent & sera not at room temperature or not well mixed before use. d) Too long time for addition of calibrator, controls, samples or reagents, Inconsistency in time intervals. E) Interference in optical pathway due to Air bubbles.	Check all 8 ports/ manifold for uniform flow of wash buffer. If there are blockage, clean the ports. Use only calibrated pipettes with well fitted tips & pipette carefully without bubbling. Equilibrate reagents to room temperature and mix thoroughly before use Develop consistent and uniform technique.
5. False Positive	a) Beside 3a, b, c, d, e & f incorrect interpretation and calculation of final results b) High incubator temperature, incorrect timing or pipetting c) Use of turbid/ lipaemic or hemolyzed sample.	Check the calculation part given in the insert and correctly interpret. Check incubator temperature, procedure & repeat assay. Centrifuge the sample at 5000 rpm for 30 minutes and re-run the test with clear sample.
6. False Negative/ low O.D. of Calibrator, Positive control/ positive sample	a) Inadequate addition of substrate/conjugate solution. b) Kit expired, reagent of different kit used. c) White particles in working substrate solution. d) Uncalibrated pipettes, improper pipetting. e) Deterioration of Enzyme conjugate, TMB Substrate/ TMB Diluent. f) Stop solution is added before 30 minutes. Reaction terminated before 30 minutes. g) O.D. taken at incorrect wavelength. h) Incorrect incubator temperature, timing or pipetting	Follow the procedure meticulously & repeat assay. Check the expiry of the kit before use. Discard the substrate and prepare the working substrate again in fresh tube. Use only calibrated pipettes with well fitted tips & pipette carefully without bubbling. Check storage of reagents. They shall be stored at 2-8°C. Follow the test procedure meticulously. Read O.D. values at 450 nm and 630 nm. Check incubator temperature, procedure & repeat assay.

