


HEPALISA

Microwell ELISA Test for the Detection of
Hepatitis B Surface Antigen (HBsAg)
in Human Serum or Plasma



- Based on "Direct Sandwich Principle"
- Detects all known 11 sub-types of hepatitis B virus
- Breakaway wells for minimal wastage
- Short Protocol: Based on Single Washing Step
- Analytical Sensitivity: 0.1 ng/ml
- Long Shelf Life: 24 months at 2-8°C
- Excellent Specificity: Increased power of discrimination between Negative & Positive samples



SENSITIVITY: 100%*
SPECIFICITY: 100%*



Scan QR Code
For more product &
service information




J. Mitra & Co. Pvt. Ltd.

.....a vision to serve mankind™


HBV: The Silent Killer

Hepatitis B is a potentially life-threatening liver infection caused by the hepatitis B virus. Originally known as "serum hepatitis."

It is a major global health problem and the most serious type of viral hepatitis which can cause chronic liver disease causing a high risk of death from cirrhosis of the liver and liver cancer. Hepatitis B virus (HBV) is a member of the Hepadnavirus family which consists of an outer lipid envelope, core and partially double stranded DNA as the genetic material. The outer surface of the virus is composed of lipid and protein, which is called the surface antigen (HBsAg), and is produced in excess during the life cycle of the virus.



Serological & Clinical Pattern: During acute HBV infection



HBsAg can be detected in high levels in serum during acute HBV or in chronic HBV whereas Anti-HBs presence usually indicates recovery & immunity from HBV infection.


IgM Anti-HBc presence indicates recent acute infection only.

The diagnosis of Hepatitis B virus (HBV) infection was revolutionized by the discovery of Australia antigen, now called Hepatitis B Surface Antigen (HBsAg).

It is the first detected during HBV infection before appearance of symptoms within 2 weeks and disappears within 3-4 months of exposure.

It is also produced in large quantities & persists for more than 6 months in carrier & chronic state with a highly antigenic determinant.

Principle



Simple Test Procedure

Add 100 μ l NC to A1 & B1, PC to C1 & D1 and sample from E1																																																																		
Prepare working conjugate solution (50x)	<table><tr><td>No. of strip</td><td>1</td><td>2</td><td>3</td><td>4</td><td>5</td><td>6</td><td>7</td><td>8</td><td>9</td><td>10</td><td>11</td><td>12</td></tr><tr><td>Enzyme conjugate</td><td>10</td><td>20</td><td>30</td><td>40</td><td>50</td><td>60</td><td>70</td><td>80</td><td>90</td><td>100</td><td>110</td><td>120</td></tr><tr><td>Concentrate (μl).</td><td>0.5</td><td>1.0</td><td>1.5</td><td>2.0</td><td>2.5</td><td>3.0</td><td>3.5</td><td>4.0</td><td>4.5</td><td>5.0</td><td>5.5</td><td>6.0</td></tr><tr><td>Conjugate Diluent(ml).</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr></table>	No. of strip	1	2	3	4	5	6	7	8	9	10	11	12	Enzyme conjugate	10	20	30	40	50	60	70	80	90	100	110	120	Concentrate (μ l).	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0	Conjugate Diluent(ml).																									
No. of strip	1	2	3	4	5	6	7	8	9	10	11	12																																																						
Enzyme conjugate	10	20	30	40	50	60	70	80	90	100	110	120																																																						
Concentrate (μ l).	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0																																																						
Conjugate Diluent(ml).																																																																		
Add 50 μ l working conjugate in each well																																																																		
Cover the plate and incubate for 60 min. at 37°C																																																																		
Wash (6 cycles)																																																																		
Prepare working substrate solution (1:1)	<table><tr><td>No. of strip</td><td>1</td><td>2</td><td>3</td><td>4</td><td>5</td><td>6</td><td>7</td><td>8</td><td>9</td><td>10</td><td>11</td><td>12</td></tr><tr><td>TMB</td><td>0.5</td><td>1.0</td><td>1.5</td><td>2.0</td><td>2.5</td><td>3.0</td><td>3.5</td><td>4.0</td><td>4.5</td><td>5.0</td><td>5.5</td><td>6.0</td></tr><tr><td>Substrate (ml).</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>TMB</td><td>0.5</td><td>1.0</td><td>1.5</td><td>2.0</td><td>2.5</td><td>3.0</td><td>3.5</td><td>4.0</td><td>4.5</td><td>5.0</td><td>5.5</td><td>6.0</td></tr><tr><td>Diluent (ml).</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr></table>	No. of strip	1	2	3	4	5	6	7	8	9	10	11	12	TMB	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0	Substrate (ml).													TMB	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0	Diluent (ml).												
No. of strip	1	2	3	4	5	6	7	8	9	10	11	12																																																						
TMB	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0																																																						
Substrate (ml).																																																																		
TMB	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0																																																						
Diluent (ml).																																																																		
Add 100 μ l working substrate in each well																																																																		
Incubate in dark for 30 min at room temperature																																																																		
Add 100 μ l Stop Solution and read result at 450 nm/630 nm																																																																		