CELIAC SELF TEST

(HOME TEST for Celiac disease Diagnosis)

Rapid Visual and sensitive one step test for the Qualitative Detection of Celiac disease associated IgM, IgG & IgA antibodies in Human Serum/Plasma

INTRODUCTION

Celiac disease (CD) is an enteropathy caused by intake of gluten proteins present in cereals such as wheat, barley and rye in genetically predisposed individuals.

Circulating IgA class anti tTG (tissue transglutaminase) and EMA (endomysium) are highly prevalent and characteristic in active celiac disease. Approximately 1% of the population develop celiac disease associated damage of gut epithelium referring to

celiac disease at varying ages, but even the majority of the cases may remain undiagnosed despite the development of chronic disease (CD) associated consequences CD can be silent, and such a condition may already lead to unfavorable sequelae.

INTENDED USE

Celiac Self Test is a visual, rapid, sensitive and accurate one step immunoassay for the qualitative detection of anti tTG IgA, IgG & IgM antibodies in Human Serum or Plasma. The assay is intended to be used as an aid in the recognition and diagnosis of Celiac disease. This test is intended to be for home use with self-collected whole blood specimens from individuals aged 18 years or above (or collected by adults for individuals below 18 years of age) as well as for professional use with collected venous Whole Blood/ Serum/ Plasma samples. This test is for in vitro diagnostic use only and is designed as an aid in the early diagnosis of celiac disease in patients with clinical symptoms.

PRINCIPLE (ANTIGEN-ANTIBODY REACTION)

Celiac Self Test is a one step immunoassay based on the antigen capture, or "sandwich" principle. The method uses tissue transglutaminase (tTG) antigen conjugated to colloidal gold and the same antigen is immobilized on a nitrocellulose strip. The diluted test sample is introduced to and flows laterally through an absorbent pad where it mixes with the signal reagent. If the sample contains celiac disease specific antibody, the colloidal gold-antigen conjugate binds to the antibody, forming an antibody-antigen-colloidal gold complex. The complex then migrates through the nitrocellulose strip by capillary action. When the complex meets the line of immobilized tTG antigen (Test line) "T", the complex is trapped forming an antigen-antibody-antigen colloidal gold complex. This forms a pink band indicating the sample is reactive for celiac disease specific antibody. Absence of a colored band in the test region indicates a non-reactive test result. A red procedural control line should always develop at 'C' region to indicate that the test has been performed properly.

KIT STORAGE AND STABILITY

1. Store the kit at 2-30ºC in cool and dry place.
2. The kit is suitable till the expiry date mentioned on the kit.
3. Do not freeze the kit.

KIT CONTENTS

Individually foil pouched test device with a dessicant. | Assay Buffer tube with nozzle | Sample Dropper, Swab & Lancet | Disposal bag | Instruction Manual

STEP-1 PREPARATION FOR THE TEST

a) Wear disposable gloves before running the test. Ensure there is no cut or wound in hand.


⚠ **DO NOT PERFORM THE TEST, in case there is a cut or wound in hand.**

b) Take out all the contents from the kit box and place it on a flat surface.

⚠ **Check the expiry date at the back of test device foil pouch . Use fresh kit if expiry date has passed.**

c) Tap the assay buffer tube gently to ensure buffer comes down to bottom of the tube. Unscrew the cap and remove from the tube. Place the nozzle on the assay buffer tube and press it gently to ensure it is tightly placed on the tube.

⚠ **Care should be taken while removing the cap so that no spillage should occur. Do not press the tube while placing the nozzle, otherwise assay buffer may come out.**

d) Take out the Celiac Self Test card from the foil pouch prior to use and place it on a flat and dry surface.

⚠ **Test should be run immediately once pouch is open.**

STEP-2 SPECIMEN COLLECTION

A) Whole blood sample collection from finger prick

a) Clean the fingertip with alcohol swab and dry completely.

b) Prick the fingertip with single use lancet provided in the kit.

⚠ **Discard the lancet in the Assay Buffer Tube once the test is completed.**

c) Collect the whole blood using sample dropper provided.

- OR -

B) Collected venous blood (Serum/ Plasma/ Whole Blood) sample from collection tubes

a)Take sample from collection tube. Fresh samples are preferred for testing as they perform best when tested immediately after collection. If samples are not immediately tested, they should be stored at 2-8°C for not more than 3 days, otherwise false / erroneous results may be obtained.


⚠ **Do not use hemolyzed/ turbid or lipaemic samples as it may lead to erroneous results.**

STEP-3 RUNNING THE TEST ON DEVICE

a) Add 4 μl of whole blood using the sample dropper into the sample well of the test card.

⚠ **Sample taken above the mark as shown in fig. (a) and sample taken below the mark as shown in fig. (c) are wrong and will lead to erratic results.**


Note: Sample dropper should be discarded immediately into the disposable bag provided after adding sample to sample well.

b) Add 2 drops assay buffer into the sample well.

⚠ **Addition of incorrect no. of drops of Assay Buffer may give invalid or false result.**

c) Allow reaction to occur for 20 minutes.


d) Read results at 20 minutes. Positive results may appear as early as 2-10 minutes. However negative results must be confirmed at 20 minutes. **Do not read results after 30 minutes.** Discard the Celiac Self test card immediately after results at 20 minutes in disposable bag considering it to be potentially infectious. Follow the local regulations for disposal of potential biohazardous material.

STEP-4 RESULT INTERPRETATION


POSITIVE RESULT: If two distinct pink line appear in the control region 'C' and test region 'T', then the specimen is reactive for Celiac antibody. Any pink line observed in test are (T) - whether faint or dark, along with Control line shall be considered as Celiac reactive.

NEGATIVE RESULT: If only one pink coloured line appear in the control region 'C', then the specimen is negative for Celiac antibody.

INVALID RESULT: If neither control line "C" nor test line "T" appears or only test line appears, the test should be treated as Invalid. Repeat the test again with new card. Invalid test may be because of following reasons:
(a) Kit is stored in high humidity and temperature.
(b) Wrong test procedure.
(c) Long atmospheric exposure of the test device after opening the pouch.

STEP-5 DISPOSAL OF USED KIT COMPONENTS

Discard the used lancet in Assay Buffer Tube; remove the nozzle from the tube and add the used lancet in tube and screw cap the tube. Discard the Assay Buffer Tube containing lancet, run device and sample dropper in disposal bag provided with the kit. Close the bag properly and throw the bag in the house-hold dustbin.

WARNING FOR USERS

1. The use of disposable gloves is RECOMMENDED while running the test.
2. Do not perform the test in area with strong air flow i.e. under fan or strong air conditioning.
3. In case there is a cut or wound in hand, DO NOT PERFORM THE TEST.
4. Do not smoke, drink or eat in areas where test is performed.
5. Mark the test card with patient's name or identification number. Improper identification may lead to wrong result reporting.
6. This test detects the presence of Celiac (tTG) antibodies in the human Whole blood/ serum/ plasma. Do not use urine and saliva as specimen for this test.
7. Dispose off all specimens and materials used to perform the test in the disposal bag provided with the kit in accordance with local regulations.
8. Wash hands thoroughly with soap or any suitable detergent after the use of the kit.

PRECAUTIONS

1. Bring kit to room temperature (20-30ºC) before use.
2. Do not use kit components beyond the expiration date, which is printed on the kit.
3. Use a separate Sample Dropper for each specimen.
4. Do not re-use the test kit.
5. Do not use the test kit if the pouch is damaged or the seal is broken.
6. Do not mix or use the assay buffer of different kits/lot.
7. Ensure that the nozzle cap is properly placed on the assay buffer tube.

LIMITATIONS AND INTERFERENCES

1. Any deviation from test procedure may lead to invalid/ erratic results.
2. This is only an in-vitro diagnostic screening test and should not be used as the sole criteria for the diagnosis of Celiac disease.
3. This is a Qualitative test and therefore neither determine Celiac (tTG) antibodies quantitative value nor the rate of Celiac (tTG) antibodies concentration.
4. A negative result may occur if the concentration of Celiac (tTG) antibodies in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly. Therefore a negative test result does not eliminate the possibility of Dengue infection, and should be confirmed by ELISA and/or Biopsy.
5. Positive test results do not rule out co-infections with other pathogens.

BIBLIOGRAPHY

1. A. George F. Davidson, Eric G. Hassal. Screening for celiac disease. Ediitorial. Can Med Assoc J 1997;157:p.547-548.
2. Jocelyn S. Fraser, Paul J. Ciclitira, Pathogenesis of celiac disease: implications for treatment, World J Gastroenterol, 2001; 7(6): 772-776.
3. T. Raivio et al, Self transglutaminase-based rapid celiac disease antibody detection by a lateral flow method, Aliment Pharmacol Ther 2006; 24: 147-154.
4. Catassi C, Rätesh IM, Fabiani E, Rossini M, Bordichhia F, Candela F, et al. Coeliac disease in the year 2000: exploring the iceberg. Lancet 1994,343:200-3
5. Korponay-Szabo IR, Dahlbom I, Laurila K, et al. Elevation of IgG antibodies against tissue transglutaminase as a diagnostic tool for coelic disease in selective IgA deficiency. Gut 2003; 52: 1567-71.

DESCRIPTION OF SYMBOLS USED

The following are graphical symbols used in or found on J. Mitra diagnostic products and packing. These symbols are the most common ones appearing on medical devices and their packing. They are explained in more detail in the European Standard EN ISO15223-1:2021.

Manufactured By	In vitro diagnostic medical device
No. of tests	Consult Instructions for use
LOT Lot Number / Batch Number	4°C 30°C Temperature Limit
Manufacturing Date	Caution, see instruction for use
Expiry Date	REF Catalogue Number
Do not use if package is damaged	Contains biological Material of Animal Origin
Single use only	Keep away from sunlight
Keep Dry	Country of Manufacture

in vitro diagnostic reagent, not for medicinal use

J. MITRA & CO. PVT. LTD.
A 180-181, Okhla Indl. Area, Phase-1, New Delhi-110 020, INDIA
Ph: +91-11-47130300, 47130500
e-mail: jmitra@jmitra.co.in Internet: www.jmitra.co.in

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