

J.MITRA & CO. PVT. LTD.

T4 QUANTI MICROLISA

QC DATA SHEET

Kit	T4 QUANTI MICROLISA		
Lot No.	ETB020724		
Mfg. Date	2024-07		
Expiry. Date	2025-12		

Standard	Concentration ( $\mu\text{g/dl}$ )	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	2.843	>1.710
S2	2.0	2.128	>1.280
S3	5.0	1.463	
S4	10.0	0.714	>0.430
S5	15.0	0.486	
S6	25.0	0.358	

Approved /  Not Approved

Quality Control

  
Signature

Date: 12-07-2024