

**J. MITRA & CO. PVT. LTD.**

**TSH QUANTI MICROLISA**

**QC DATA SHEET**

Kit	<b>TSH QUANTI MICROLISA</b>
Lot No.	ETH040724
Mfg. Date	2024-07
Expiry Date	2026-06

Standard	Concentration ( $\mu$ IU/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	0.013	
S2	0.5	0.099	
S3	2.5	0.383	>0.230
S4	5.0	0.667	
S5	10.0	1.370	>0.820
S6	20.0	2.208	
S7	40.0	2.725	>1.640

✓  
Approved / ~~Not Approved~~

Quality Control

  
Signature

Date 30-07-2024