

**J.MITRA & CO. PVT. LTD.**


**PCT QUANTI MICROLISA**

**QC DATA SHEET**

Kit	<b>PCT QUANTI MICROLISA</b>
Lot No.	EPT010824
Mfg. Date	2024-08
Expiry. Date	2026-01

Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range
S1	0.0	0.010	
S2	0.5	0.066	>0.020
S3	1.0	0.185	
S4	2.5	0.437	>0.340
S5	10.0	1.294	
S6	25.0	2.458	>1.630

✓  
Approved / Not Approved

  
Signature

Date 05-08-2024