



**J. Mitra & Co. Pvt. Ltd.**

*.....a vision to serve mankind™*

## ELISA WASHER



## ELISA READER



### ELISA WASHER

#### ***Features:***

- Provided with 3 wash bottles for switching off to 3 different ELISA tests.
- In built Incubator adjustable to Required temperature(Ambient, 30°C& 37°C).
- Provision for 2 ELISA PLATES (96T) incubation at a time
- Inbuilt shaker provided for better washing
- Accurate results due to cover lid as no Air Draft at the time of sample and Reagents incubation.
- Provision of auto rinsing of tubes/probes & hence No cross contamination

### ELISA READER

#### ***Features:***

- User friendly to operate with Touch screen display.
- Easily programmable for Qualitative and Quantitative tests.
- Includes 4PL Method for Quantitative tests.
- In-built thermal printer for printing test results
- Storage Capacity: 1,000 patient data or 10,000 sample records.

# WIDE RANGE OF QUANTITATIVE/ QUALITATIVE ELISA TESTS

## .....23 tests to choose from!

- Excellent Sensitivity & Specificity for all qualitative tests
- High precision, good reproducibility, linearity & accuracy co-efficient ( $R^2$ ) of > 0.98 for all quantitative tests

DENGUE NS1Ag MICROLISA



DENGUE IgM MICROLISA



DENGUE IgG MICROLISA



CHIKUNGUNYA IgM MICROLISA



LEPTO IgM MICROLISA



MICROLISA HIV Ag & Ab



MICROLISA HIV



HIV MICROLISA 3.0



HCV MICROLISA



TSH QUANTI MICROLISA



T4 QUANTI MICROLISA



T3 QUANTI MICROLISA




MALARIA Ag MICROLISA



SCRUB TYPHUS IgM MICROLISA



HEPALISA ULTRA



HEPALISA



COVID KAWACH IgG MICROLISA



COVID 19 (IgM+IgG+IgA) Microlisa



CELIAC MICROLISA



hs-CRP QUANTI MICROLISA



PCT QUANTI MICROLISA



VITAMIN B12 QUANTI MICROLISA



VITAMIN D QUANTI MICROLISA



Also Connect with us on:



ISO 9001:2015



ISO 13485:2016  
ICMED 13485



For further enquiries, Please contact:  
**J. Mitra & Co. Pvt. Ltd.**

A 180-181, Okhla Indl. Area, Phase-1, New Delhi - 110 020, INDIA  
Tel: +91-11-47130300, 47130500  
E-mail: [jmitra@jmitra.co.in](mailto:jmitra@jmitra.co.in) Website: [www.jmitra.co.in](http://www.jmitra.co.in)