



**J. Mitra & Co. Pvt. Ltd.**

*.....a vision to serve mankind™*

J. MITRA OFFERS


STATE OF THE ART TECHNOLOGY

# **ELISA READER & ELISA WASHER**

*28 tests to choose from.....*



# Offering Excellent Quality & Trust through Wide Range of Quantitative & Qualitative ELISA Tests...



# ER-181s Microplate Reader



In-built thermal printer for printing test results



0 - 4.00 Abs  
Enhanced Measuring Range



Touch screen & 5.7" LCD display



Storage Capacity: 100 test protocols (28 pre-programmed protocols) and 10,000 Sample records



8-channel optical fiber system (more accurate results)



Multi-assay: Up to 12 different assays on one plate



Plug-n-play usage with PC



Lamp sleeping and wake-up function (cost saving & increase life of the lamp)



HIS & LIS compatible with application support



Auto self-check at start-up



Data Protection (Auto back up in case of power outage)



Online technical support for programming/trouble shooting



User friendly Windows Operating System



# EW-181s

## Microplate Washer

**NEW**



Built-in incubator with user adjustable temperature settings at 25°C, 30°C and 37°C



Provision for incubation of two 96 well plates at a time and alarm system



Multi-channel, 3 bottles: 1 each for wash, rinse & waste



5" LCD display



Large memory, upto 100 wash protocols can be programmed




Plate shaking function, time and speed are adjustable



Online technical support for programming/trouble shooting



Liquid level detection and alert function



8-way manifold, automatic & manual positioning



User friendly operating system



Wash cycle monitoring



Easy maintenance

Also Connect with us on:



ISO 13485:2016  
ICMED 13485



For further enquiries, Please contact:

**J. Mitra & Co. Pvt. Ltd.**

A 180-181, Okhla Indl. Area, Phase-1, New Delhi - 110 020, INDIA

Tel: +91-11-47130300, 47130500

E-mail: [jmitra@jmitra.co.in](mailto:jmitra@jmitra.co.in) Website: [www.jmitra.co.in](http://www.jmitra.co.in)