

J.MITRA & CO. PVT. LTD.


T3 QUANTI MICROLISA

QC DATA SHEET

Kit	T3 QUANTI MICROLISA
Lot No.	ETA020325
Mfg. Date	2025-03
Expiry. Date	2026-08

Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	2.405	>1.440
S2	0.5	2.055	
S3	1.0	1.576	>0.950
S4	2.5	1.024	
S5	5.0	0.674	>0.400
S6	7.5	0.433	

✓
Approved / ~~Not Approved~~
Quality Control


Signature
Date : 26-03-2025