

**J. MITRA & CO. PVT. LTD.**

**TSH QUANTI MICROLISA**

**QC DATA SHEET**

|             |                             |
|-------------|-----------------------------|
| Kit         | <b>TSH QUANTI MICROLISA</b> |
| Lot No.     | ETH040925                   |
| Mfg. Date   | 2025-09                     |
| Expiry Date | 2027-08                     |

| Standard | Concentration<br>( $\mu$ IU/ml) | Avg. O.D. at<br>450/630 nm | Acceptable range<br>of O.D. |
|----------|---------------------------------|----------------------------|-----------------------------|
| S1       | 0                               | 0.018                      |                             |
| S2       | 0.5                             | 0.083                      |                             |
| S3       | 2.5                             | 0.287                      | >0.170                      |
| S4       | 5.0                             | 0.668                      |                             |
| S5       | 10.0                            | 1.285                      | >0.770                      |
| S6       | 20.0                            | 1.933                      |                             |
| S7       | 40.0                            | 2.600                      | >1.560                      |

Approved / Not Approved

Quality Control

  
Signature

Date: 10-09-2025