

**J. MITRA & CO. PVT. LTD.**

**TSH QUANTI MICROLISA**


**QC DATA SHEET**

Kit	<b>TSH QUANTI MICROLISA</b>
Lot No.	ETH010326
Mfg. Date	2026-03
Expiry Date	2028-02

Standard	Concentration ( $\mu$ IU/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	0.017	
S2	0.5	0.086	
S3	2.5	0.308	>0.180
S4	5.0	0.721	
S5	10.0	1.189	>0.710
S6	20.0	2.031	
S7	40.0	2.586	>1.550

✓  
Approved / Not Approved-

Quality Control

  
Signature

Date: 25-03-2026