

**J. MITRA & CO. PVT. LTD.**

**VITAMIN D QUANTI MICROLISA**

**QC DATA SHEET**

Kit	VITAMIN D QUANTI MICROLISA
Lot No.	EVD010526
Mfg. Date	2026-05
Expiry Date	2027-10


Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range
S1	0	2.487	>1.490
S2	4	2.189	
S3	10	1.911	>1.150
S4	25	1.104	
S5	60	0.597	
S6	130	0.098	

Standard	Concentration (ng/ml)	Acceptable range (ng/ml)	Observed Concentration (ng/ml)
Low	20	15-25	17.49
High	50	40-60	49.95

✓  
Approved / ~~Not Approved~~

Quality Control

Signature



Date:- 06-05-2026