

Kit	VITAMIN B12 QUANTI MICROLISA
Lot No.	EVB010726
Mfg. Date	2026-07
Expiry Date	2027-12

Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	2.147	>1.290
S2	100	1.749	
S3	200	1.460	>0.880
S4	400	1.113	
S5	1000	0.590	>0.350
S6	2000	0.367	

Approved / Not Approved

Quality Control



Signature

Date: 24-07-2026