

Kit	TSH QUANTI MICROLISA
Lot No.	ETH020626
Mfg. Date	2026-06
Expiry Date	2028-05

Standard	Concentration (μ IU/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	0.017	
S2	0.5	0.089	
S3	2.5	0.331	>0.200
S4	5.0	0.717	
S5	10.0	1.189	>0.710
S6	20.0	2.064	
S7	40.0	2.579	>1.550

Approved / Not Approved

Quality Control



Signature

Date: 16-06-2026