

T3 QUANTI MICROLISA

QC DATA SHEET

Kit	T3 QUANTI MICROLISA
Lot No.	ETA030419
Mfg. Date	2019-04
Expiry. Date	2020-09

Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	2.534	>1.567
S2	0.5	2.012	
S3	1.0	1.734	>1.104
S4	2.5	1.089	
S5	5.0	0.643	>0.467
S6	7.5	0.357	

QC Approved / Not Approved

Quality Control


Signature

Date 12/04/2019