

J.MITRA & CO. PVT. LTD.

T4 QUANTI MICROLISA

QC DATA SHEET

Kit	T4 QUANTI MICROLISA
Lot No.	ETB010322
Mfg. Date	2022-03
Expiry. Date	2023-08

Standard	Concentration ($\mu\text{g}/\text{dl}$)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	2.710	>1.630
S2	2.0	2.172	>1.300
S3	5.0	1.322	
S4	10.0	0.657	>0.390
S5	15.0	0.458	
S6	25.0	0.246	

 ~~Approved / Not Approved~~

Quality Control

Signature



Date 25-03-2022