

J.MITRA & CO. PVT. LTD.

TSH QUANTI MICROLISA

QC DATA SHEET

Kit	TSH QUANTI MICROLISA
Lot No.	ETH040622
Mfg. Date	2022-06
Expiry. Date	2024-05

Standard	Concentration (μ IU/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	0.012	
S2	0.5	0.075	
S3	2.5	0.230	>0.140
S4	5.0	0.526	
S5	10.0	1.194	>0.720
S6	20.0	2.196	
S7	40.0	2.737	>1.640

Approved / Not Approved

Quality Control


Signature

Date 10/06/2022