

J.MITRA & CO. PVT. LTD.

VITAMIN D QUANTI MICROLISA

QC DATA SHEET

Kit	VITAMIN D QUANTI MICROLISA
Lot No.	EVD031023
Mfg. Date	2023-10
Expiry. Date	2025-03

Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range
S1	0	2.358	>1.410
S2	4	1.879	
S3	10	1.643	>0.990
S4	25	1.225	
S5	60	0.591	
S6	130	0.105	

Standard	Concentration (ng/ml)	Acceptable range (ng/ml)	Observed Concentration(ng/ml)
Low	20	15-25	21.56
High	50	40-60	53.43

Approved / Not Approved

Quality control


Signature

Date 01-11-2023