

J.MITRA & CO. PVT. LTD.

VITAMIN D QUANTI MICROLISA

QC DATA SHEET

Kit	VITAMIN D QUANTI MICROLISA
Lot No.	EVD030824
Mfg. Date	2024-08
Expiry. Date	2026-01

Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range
S1	0	2.717	>1.630
S2	4	2.114	
S3	10	1.737	>1.040
S4	25	1.237	
S5	60	0.597	
S6	130	0.089	

Standard	Concentration (ng/ml)	Acceptable range (ng/ml)	Observed Concentration(ng/ml)
Low	20	15-25	20.50
High	50	40-60	53.10

✓
Approved / ~~Not Approved~~

Quality control

Signature 

Date 21-08-2024