

J. MITRA & CO. PVT. LTD.

TSH QUANTI MICROLISA


QC DATA SHEET

Kit	TSH QUANTI MICROLISA
Lot No.	ETH050924
Mfg. Date	2024-09
Expiry Date	2026-08

Standard	Concentration (μ IU/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	0.014	
S2	0.5	0.086	
S3	2.5	0.305	>0.180
S4	5.0	0.683	
S5	10.0	1.308	>0.780
S6	20.0	2.080	
S7	40.0	2.753	>1.650

✓
Approved / ~~Not Approved~~

Quality Control


Signature

Date 28-09-2024