

J.MITRA & CO. PVT. LTD.


T3 QUANTI MICROLISA

QC DATA SHEET

Kit	T3 QUANTI MICROLISA
Lot No.	ETA041224
Mfg. Date	2024-12
Expiry. Date	2026-05

Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	2.522	>1.510
S2	0.5	2.021	
S3	1.0	1.560	>0.940
S4	2.5	0.957	
S5	5.0	0.592	>0.350
S6	7.5	0.290	

✓
Approved / ~~Not Approved~~
Quality Control


Signature

Date : 06-12-2024