

J.MITRA & CO. PVT. LTD.

T3 QUANTI MICROLISA

QC DATA SHEET

Kit	T3 QUANTI MICROLISA
Lot No.	ETA010125
Mfg. Date	2025-01
Expiry. Date	2026-06

Standard	Concentration (ng/ml)	Avg. O.D. at 450/630 nm	Acceptable range of O.D.
S1	0	2.561	>1.540
S2	0.5	2.164	
S3	1.0	1.600	>0.960
S4	2.5	0.923	
S5	5.0	0.534	>0.320
S6	7.5	0.310	

✓
Approved / ~~Not Approved~~
Quality Control

Signature

Date : 08-01-2025